

OXYPLETH[®]

Service Manual

Pulse Oximeter

Model 520A

June 5, 2001

Catalog Number 5693-90-01

Novamatrix Medical Systems Inc. Wallingford, Connecticut, U.S.A. 06492.
Copyright ©1992-2001. All rights reserved. No part of this manual may be reproduced
without the written permission of Novamatrix Medical Systems Inc.

Guarantee

Equipment manufactured or distributed by Novamatrix Medical Systems Inc., is fully guaranteed, covering materials and workmanship, for a period of one year from the date of shipment, except for certain disposable products and products with stated guarantees other than one year. Novamatrix reserves the right to perform guarantee service(s) at its factory, at an authorized repair station, or at the customer's installation.

Novamatrix' obligations under this guarantee are limited to repairs, or at Novamatrix' option, replacement of any defective parts of our equipment, except fuses, batteries, and calibration gasses, without charge, if said defects occur during normal service.

Claims for damages during shipment must be filed promptly with the transportation company. All correspondence concerning the equipment must specify both the model name and number, and the serial number as it appears on the equipment.

Improper use, mishandling, tampering with, or operation of the equipment without following specific operating instructions will void this guarantee and release Novamatrix from any further guarantee obligations.

Service Department
For factory repair service, call toll free
1-800-243-3444
In Connecticut, call Collect (203) 265-7701
Facsimile (203) 284-0753
World Wide Web: <http://www.novamatrix.com>
Internet: techline@novamatrix.com

Caution: Federal (U.S.A.) law restricts this device to sale, distribution, or use by or on the order of a licensed medical practitioner.

Copyright© 1992-2001 Novamatrix Medical Systems Inc. This document contains information which is proprietary and the property of Novamatrix Medical Systems Inc., and may not be reproduced, stored in a retrieval system, translated, transcribed, or transmitted, in any form, or by any means, without prior explicit written permission from Novamatrix Medical Systems Inc.

Novamatrix and OXYPLETH are a registered trademarks. SuperBright is a trademark of Novamatrix Medical Systems Inc.

The OXYPLETH® monitor and its sensors and accessories are covered by the following US patents: 5,190,038 5,398,680 5,448,991 5,820,550 5,999,834 5,891,026 6,073,038 6,149,481. Other patents pending

Revision History

August 31, 1992 Release 00. This manual replaces all preliminary versions.
Based on Revision 2.0 CPU software.

June 5, 2001 Revision 01. R-N905


Declaration of Conformity with European Union Directives

The authorized representative for Novamatrix Equipment is:

European Compliance Services Limited
Oakdene House
Oak Road
Watchfield
Swindon, Wilts SN6 8TD
UK

Manufacturing, Quality and Safety

Novamatrix manufacturing facility is certified to ISO 9001 and EN46001 (MDD93/42/EEC Annex II). Novamatrix Medical Systems Inc. products bear the “CE 0086” mark. The product is certified by Underwriter’s Laboratories (UL) to bear the UL mark; and tested by TUV Rheinland to IEC601-1 / EN60601-1.



Contents

Contents	iii
List of Figures	vii
List of Tables.....	viii
Introduction.....	1
Purpose	1
Technology Description	1
Conventions Used In This Manual.....	2
Acknowledgments	2
Patient Safety	3
Warnings.....	4
Cautions.....	5
Front Panel	6
Rear & Top Panels	7
Summary of Operation	8
Menu Trees.....	9
Electronic Theory of Operation	15
2472 Power Supply Board.....	15
AC Mains and Battery Operation Overview	15
AC Mains Operation	16
Battery Operation.....	16
2726 Power Supply Board.....	17
AC Mains and Battery Operation Overview	17
AC Mains Operation	17
Battery Operation.....	18
2471 & 2775 Main Board.....	19
Power On/Off Control Circuitry	19
Power Supplies.....	20
Voltage References	20
Preserving RAM and Real Time Clock Data.....	21

Low Battery Voltage Shutdown.....	21
Timing Sequencer.....	22
Data Sampling Controller.....	22
Sensor LED Drive Circuits	23
Sensor Photodiode Return Path	24
Calibrating the 20-Bit Analog-to-Digital Convertors	25
20-Bit Analog-to-Digital Conversion.....	26
Sensor Status Decoding and Conversion.....	26
Sensor Status Parameters.....	27
Microprocessor	28
Memory	28
Real Time Clock (RTC).....	29
Sound generator	29
Keypanel Interface.....	30
Display Interface	30
I/O Device Controller	31
Watchdog Timer	31
Serial I/O Controller	32
RS232 Serial Communication.....	33
Maintenance	35
General.....	35
Maintenance Schedules	35
Cleaning and Sterilization.....	36
Model 520A Monitor.....	36
Finger Sensor	36
Y-SENSOR™ and Y-STRIP™ Taping System.....	36
Battery Life and Maintenance	37
Mains Voltage Configuration	37
Fuse Replacement.....	37
Changing the Mains Voltage Setting	38
AC Mains on older style assemblies.....	39
Replacing the AC Mains Fuse(s)	40
Changing the AC Mains Voltage.....	40
Assembly Exchanges	41
Changing System Software	44
Changing System Software on 2471 Main Board.....	44
Troubleshooting.....	47
Functional Test	51
Introduction	51
Monitor Functional Test	51
Special Power Up Functions	53

Accuracy Test	55
Introduction	55
Monitor Accuracy Test	55
Calibration Tests.....	59
Introduction	59
Equipment Required and Test Setup	59
Test Procedure	60
Safety Testing.....	63
Connecting to other Equipment	65
Connecting the Seiko DPU-414 Thermal Printer	65
Configuring the Seiko DPU-414 Printer	65
Setting the DIP switches:.....	65
Connecting the ThinkJet Printer	67
Connecting the Model 315 Printer	67
Connecting Seiko DPU-411 Thermal Printer	68
Analog Output Module	69
Analog Output Setup	70
Specifications.....	73
General	73
Oxygen Saturation (SpO2) Section	73
Pulse Rate Section	73
General Specifications.....	74
Additional Features.....	74
Accessories.....	76
Model 520A	76
Parts Lists.....	81
Family Tree.....	81
5693-01 MAIN ASSEMBLY	81
6538-01 REAR PANEL ASSY	82
5673-01 REAR PANEL ASSY	82
2726-01 POWER SUPPLY BD ASSY	83
2472- 01 POWER SUPPLY BOARD	84
5719-01 FRONT PANEL ASSY	85
2473-01 ALERT BOARD ASSY	85
5720- 01 DISPLAY ASSY.....	85
5713-01 SPEAKER ASSY	85
5714-01 BATTERY HARNESS ASSY	85
5728-01 CABLE ASSY	86
5765-01 TOP COVER ASSY	86
5766-01 BOTTOM COVER ASSY	86

2755-01 TONE GEN REPL BOARD	87
2775-01 MAIN BOARD ASSY	87
2471-28 MAIN BD KIT (part of 2775-01)	89
2471-01 MAIN BOARD ASSY	91
Schematic and Assembly Drawings	93
Notes	95

F

List of Figures

Figure 1. Front Panel Description	6
Figure 2. Rear and Top Panel Descriptions	7
Figure 3. Figure 1. Front-End Timing	23
Figure 4. Power Entry Fuse Access Door Opening	40
Figure 5. Fuse Removal	40
Figure 6. Power Entry Module Voltage Selector Drum Removal	41
Figure 7. Power Entry Module Voltage Selection Adjustment	41
Figure 8. Removing the Top Cover	42
Figure 9. Monitor Assembly	42
Figure 10. 520A Assembly with 2471 Main Board Removed	43
Figure 11. Changing the System Software EPROM	45
Figure 12. Softkey Identification	54
Figure 13. Connecting the ThinkJet Printer	67
Figure 14. Connecting the Model 315 Printer	68
Figure 15. Configuration of Seiko Thermal Printer	68
Figure 16. Connecting the Seiko Model DPU-411 Printer.	69
Figure 18. Analog Output Module Pinouts	70
Figure 17. Analog Output Module	70
Figure 19. Analog Output Setup	71



List of Tables

Table 1.	Display Command/Data table.....	30
Table 2.	CSI/O Decode Lines.....	32
Table 3.	SpO2 Display tolerances for TB500B settings.....	57
Table 4.	Analog Output Module Voltage Tolerances.....	71
Table 5.	Schematic and Assembly Drawings.....	93

1

Introduction

Purpose

1.1

This manual has been prepared for technicians servicing the Novametrix Model 520A Pulse Oximeter monitor. It presents technical information relating to the monitor's theory of operation, maintenance, calibration and repair. Refer to the Pulse Oximeter Model 520A User's Manual (Catalog Number 5693-23) for additional information.

Technology Description

1.2

Pulse oximetry is a non-invasive method of monitoring the oxygen saturation of arterial blood. Oxygen saturation monitoring is intended to be used in a variety of clinical situations including, respiratory therapy, anesthesia, the intensive care unit (ICU) and neonatal (NICU) and pediatric (PICU) intensive care units.

The Model 520A Pulse Oximeter provides reliable continuous measurement, display, and alerts for oxygen saturation (SpO₂) and pulse rate. The monitor can be powered from its rechargeable internal battery or from the AC Mains.

The Model 520A utilizes sensors containing two light emitting diodes (LEDs) and a photodiode. Each LED emits a specific wavelength of light (660 and 940 nanometers) through a pulsating vascular bed to the photodiode. Oxygen saturated blood absorbs different amounts of light at each wavelength as compared to unsaturated blood. Therefore, the amount of light absorbed by the blood in each pulse can be used to calculate saturation.

The Model 520A is calibrated to display "functional" saturation. This differs from the "fractional" saturation value displayed by most co-oximeters. Functional saturation is defined as:

$$\text{Functional Saturation} = \frac{\text{HbO}_2}{100 - (\text{COHb} + \text{METHb})}$$

HbO₂ = Fractional Hemoglobin

COHb = Carboxyhemoglobin

METHb = Methemoglobin

This can be considered to represent the amount of oxyhemoglobin as a percentage of the hemoglobin that can be oxygenated. Dysfunctional hemoglobins (COHb and METHb) are not included.

Pulse Rate is calculated by measuring the time interval between detected peaks of the infrared light waveform. The inverse of this measurement is displayed as pulse rate.

The oxygen saturation and pulse rate values are displayed on monitor's displays. The displayed values are updated once each second. Presence of a pulse is shown as a waveform on the display and indicated audibly by a user selectable "beep".

The Model 520A must be used in conjunction with SuperBright™ Sensors. These sensors have an 8700 series part number (e.g., 8776 or 8791). An INCOMPATIBLE PROBE display message will indicate a non-SuperBright™ Sensor (e.g., 86xx series) is in use.

Conventions Used In This Manual

1.3

The following conventions will be used throughout this manual:

- Normal text will be shown in this type.
- Message Center alerts and displays will be shown in this type.
- The names of the front panel pushbuttons (keys) will be shown **in this type**.

Acknowledgments

1.4

SuperBright, Y-STRIP and Y-SENSOR are trademarks of Novamatrix Medical Systems Inc. SARAcap is a registered trademark of Allegheny International Medical Technology, Inc. (PPG Biomedical Systems).

2

Patient Safety

The *OXYPLETH* Pulse Oximeter Monitor, Model 520A, SpO₂ Input is electrically isolated. Patient leakage current flowing from the instrument to ground is limited to less than 10 μ A at 120 VAC, 60 Hz. Patient isolation is greater than 10 M Ω , 2500 VAC rms at 60 Hz.

For maximum patient and operator safety, the following are recommended:

- **Failure of Operation:** If the monitor fails to respond as described, do not use it until the situation has been corrected by qualified personnel.
- Keep *OXYPLETH* and its accessories clean.
- Do not operate *OXYPLETH* when it is wet due to spills or condensation.
- Do not operate *OXYPLETH* if it appears to have been dropped or damaged.
- Connect the line cord only to a grounded hospital-grade outlet. *OXYPLETH* should be connected to the same electrical circuit as other equipment in use on the patient. Outlets on the same circuit can be identified by the hospital's engineering department.
- Care should be exercised to assure continued peripheral perfusion distal to the SpO₂ sensor site after application.
- Components of this product and its associated accessories which may have patient contact are free of latex.

3

Warnings



WARNING

Indicates a potentially harmful condition that can lead to personal injury

- **Explosion Hazard:** Do NOT use *OXYPLETH* in the presence of flammable anesthetics. Use of this instrument in such an environment may present an explosion hazard.
- **Electrical Shock Hazard:** Always turn *OXYPLETH* off and remove line cord before cleaning it. Do NOT use a damaged sensor or one with exposed electrical contacts. Refer servicing to qualified service personnel.
- Do not operate *OXYPLETH* when it is wet due to spills or condensation.
- Do not operate *OXYPLETH* if it appears to have been dropped or damaged.
- **Patient Safety:** Extreme care should be exercised with neonates to assure continued circulation distal to the sensor site after application.
- **Failure of Operation:** If the monitor fails to respond as described, do not use it until the situation has been corrected by qualified personnel.
- **Patient Safety:** Care should be exercised to assure continued peripheral perfusion distal to the SpO₂ sensor site after application.
- **Data Validity:** Inaccurate SpO₂ and/or Pulse Rate measurements can be caused by any of the following:
 - Incorrect application or use of a sensor
 - Significant levels of dysfunctional hemoglobin such as carboxyhemoglobin or methemoglobin
 - Significant levels of indocyanine green, methylene blue, or other intravascular dyes
 - Exposure to excessive illumination such as surgical lamps—especially ones with a xenon light source, or direct sunlight
 - Excessive patient movement, venous pulsations, electrosurgical interference
- **Data Validity:** The Pulse Oximeter should not be used as a substitute for an ECG monitor. The oximeter's Pulse Rate display reflects the pulsatile flow found at the patient extremity connected to the sensor. This rate can be affected by many factors and may occasionally be “frozen.”
- **Do NOT** attach an SpO₂ sensor distal to a blood pressure cuff. Valid data **CANNOT** be processed when the cuff is inflated. Attach the sensor to the limb opposite to the site used for the blood pressure cuff.
- **Do NOT** apply Y-Sensor tapes or wraps so tightly that the circulation is restricted. Inspect site often for adequate circulation - at least once every four hours. When applying sensors take note of the patient's physiological condition. For example, burn patients may exhibit more sensitivity to heat and pressure and therefore additional consideration such as more frequent site checks may be appropriate.

4

Cautions


CAUTION

Indicates a condition that may lead to equipment damage or malfunction.

- Do not operate *OXYPLETH* when it is wet due to spills or condensation.
- Do not operate *OXYPLETH* if it appears to have been dropped or damaged.
- Never sterilize or immerse the monitor in liquids.
- Do not sterilize or immerse sensors except as directed in this manual.
- Tension should not be applied to the sensor cable.
- Overstretching the pulse oximeter finger sensor can damage the sensor and potentially affect pulse oximeter readings. Do not stretch the finger sensor open beyond the limit for which it was designed. Overstretching can be prevented: avoid opening the sensor by any means other than squeezing the grips; Do *NOT* force the sensor onto large objects such as the bed rail.
- Do not store the monitor or sensors at temperatures less than 14° F (-10° C) or greater than 131° F (55° C).
- Do not operate the monitor or sensors at temperatures less than 50° F (10° C) or greater than 104° F (40° C).
- Federal (U.S.A.) law restricts this device to sale, distribution, or use by or on the order of a licensed medical practitioner.

5

Front Panel




1. **Saturation% Display Area;** SpO₂ and alert limit settings are displayed here.
2. **Pulse Rate Display Area;** Pulse Rate and alert limit settings are displayed here.
3. **Message Center;** area where system messages are displayed. The functions of the softkeys (#4) are annotated here.
4. **Softkeys # 1-5;** Softkeys 1-5 (left to right) cause the action annotated in the lower half of the Message Center to occur.
5. **Carrying Handle;** monitor carrying handle molded into case
6. **Sensor Input Connector;** Connect SuperBright™ Sensors here. **Press tab on sensor connector to remove sensor. Do not twist connector.**
7. **Two Minute Silence Indicator;** Illuminates (yellow) when the **AUDIO** key is pressed. SpO₂ and Pulse Rate alarms are silenced for two minutes.
8. **Audio key;** Press and release **AUDIO** to turn on/off the two minute silence function. Press and hold (approx. 3 seconds) **AUDIO** to enable the Audio Off feature (unless disabled via Options Menu). Press and release to disable Audio Off.
9. **Audio Off Indicator;** Flashes (yellow) as a warning that the audible alarms have been disabled.
10. **Alert Indicator;** Flashes (red) when an alert/alarm occurs. Continues to flash until condition corrected and **ALERT RESET** is pressed.
11. **Alert Reset key;** Press **ALERT RESET** to disable any active alert indicators. Alerts will reactivate if alert condition still exists.
12. **AC Power Indicator;** Illuminates (green) when the monitor is connected to an AC (Mains) power source and the rear panel power switch is set to "I".
13. **Power key;** Press **POWER** to turn the monitor off and on.
14. **Low Battery Indicator;** Illuminates (red) when the monitor is powered from its internal battery and less than 30 minutes of battery power remain.
15. **Red Alert Bar;** Flashes (red) when an alert/alarm occurs (unless disabled via Options Menu). Continues to flash until condition corrected and **ALERT RESET** is pressed (unless "unlatched" by the user via the Options Menu).
16. **Front Feet;** Rubber tipped front feet (2). **Kickstand;** Two position kickstand (not shown) lifts front of monitor for viewing from above.
17. **Contrast;** Press and hold for display contrast adjustment, release when desired contrast is obtained.
18. **Event;** Press to mark an event in trend memory.

Figure 1. Front Panel Description

6

Rear & Top Panels



1. **Line Cord Clip:** This clip can be set around the line cord strain relief so that the cord cannot be pulled out of the connector.
2. **Line Cord Connector:** The AC (Mains) line cord attaches to the monitor here.
3. **Fuse Compartment:** The AC (Mains) line fuse(s) are inside this compartment. Pry open with small screwdriver.
4. **AC Mains Voltage:** The currently selected AC Mains input voltage is identified here.
5. **AC Mains Power Switch:** With switch in "O" position, AC Mains voltage does not enter monitor. With switch in "I" position, AC Mains voltage allowed into monitor to power unit and/or charge internal battery.
6. **Rear Feet:** Rubber tipped rear feet (2).
7. **Serial Output Connector:** Serial (RS232) data output here for use with Tele-Sat™ telemetry system, optional analog output module, and other RS232 interfaces. A female 25-pin "D" connector serves as the interface connector.
8. **Top Cover**
9. **Carrying Handle:** Monitor carrying handle molded into case.
10. **Warning Label:** Explosion and electrical shock warnings.
11. **Patient Isolation Label:** The Model 520A is Type BF equipment.

Figure 2. Rear and Top Panel Descriptions

7

Summary of Operation


This section summarizes the operation of the Model 520A Pulse Oximeter. It is intended as a quick reference and refresher for persons who have thoroughly reviewed the Model 520A User's Manual (PN: 5693-23). Persons unfamiliar with the Model 520A should thoroughly examine the User's Manual before referencing the steps listed here.


1. If powering the monitor from the AC line voltage; Connect the line cord to the monitor and plug into a properly grounded 3-wire outlet. Set the rear panel power (Mains) switch to the I (on) position.
2. Press the front panel **POWER** key to turn the monitor on. Verify the Message Center displays Connect SpO2 Probe after the self-tests are complete.
3. Use the Menu System to modify the operation of the monitor to best suit your application. The Menu System menus include:
 - **ALRT** softkey - Use **SEL** to select SpO2 or Pulse limits, then use the ↑ or ↓ to set the limit to the desired value. Note that the monitor will maintain a spread of at least 5 digits between the upper and lower limit values.
 - **TRND** softkey - enters Trend mode and draws trend memory on display.
 - **MENU** then **AUDIO** softkeys - Select **Pulse** to turn on/off and control the volume of the audible beep with each detected pulse. Select **ALERT** to control alert volume.
 - **MENU** then **AVG** - Select either 2 second or 8 second SpO2 averaging time.
 - **MENU** then **LITE** softkeys - toggles backlight between bright and dim.
4. Verify that the SpO₂ and Pulse Rate alert limit settings are consistent with your monitoring application.
5. Connect a SuperBright™ (87xx series) sensor to the front panel input. Apply the sensor to the patient.
6. Check that within several seconds the patient's SpO₂ and Pulse Rate are displayed on the monitor.
7. As necessary, refer to the appropriate sections within the User's Manual for detailed operating instructions and explanations.


8

Menu Trees

The *OXYPLETH* menus are described on the following pages.








* The NOVACARD option must be programmed in MONITOR OPTIONS 2 menu.

A serial printer interface must be selected in MONITOR OPTIONS 2 menu.

WAVEFORM AUTO SIZE must be turned ON in the SPO2 SETUP OPTIONS menu.

* PRINT and SIZE will appear only if selected in the setup menu. CARD will appear if NOVACARD interface is selected.

1



9

Electronic Theory of Operation

The electronic theory of operation of the Model 520A Pulse Oximeter monitor is detailed in the subsections below. Refer to Section 19, *Schematic and Assembly Drawings*, on page 93 for further information.

2472 Power Supply Board

9.1

The 2472 Power Supply Board contains the circuitry needed to power the monitor from the AC Mains (line voltage). The power supply board also connects to the battery and contains the battery charging circuitry.

AC Mains and Battery Operation Overview

9.1.1

The Model 520A can be powered from its internal 12 volt battery or from the AC Mains. The green \sim (sine wave shaped) front panel indicator illuminates when the line cord is connected and the rear panel power switch is in the “I” (On) position. This indicates that AC Mains power is reaching the power supply, that the battery is being charged, and that if the monitor is turned on, it is being powered from the line voltage.

If AC Mains power is removed by unplugging the line cord or setting the rear panel power switch to the “0” (Off) position, the monitor will operate for up to four hours from its internal 12 volt lead-acid battery. As the battery voltage runs low (≈ 11.5 volts), the red E (battery indicator) on the front panel illuminates. At this point, the AC Mains should be reconnected to power the monitor and charge the battery.

If the monitor continues to be powered from a battery in a low voltage situation, at approximately 11 volts, a continuous alarm sounds for thirty seconds while the Message Center displays BATTERY EXHAUSTED CONNECT LINE CORD. If this alarm/message is ignored, the monitor displays will shut down and the battery indicator will flash on and off about every 5 seconds. If AC power is now restored, the monitor will re-initialize (go through the power up and self-test routines) and resume normal operation. However, continued battery operation will eventually activate a hardware low voltage circuit (≈ 8.5 volts) that shuts the monitor off to prevent damage to the battery. Once the unit is shut down with the hardware circuit,

the AC Mains must be connected and the front panel **POWER** key pressed before the monitor will turn back on.

AC Mains Operation

9.1.2

The AC Mains voltage enters the monitor at the rear panel Power Entry Module (PEM). This device contains a built in RFI power line filter, a double-pole single-throw switch that opens and closes both AC input lines, fuses, and an input voltage selection card.

The filtered, switched and fused output of the Power Entry Module is fed to the primary coils of the rear panel mounted system transformer, T301. The secondary output from T301 is rectified by D1 (bridge rectifier) and filtered by C1. The (loaded) DC voltage at this point is approximately 20 volts.

The 20 DC volts is fed to the 2471 main board through fuse F301 to connector E302, and is switched to the battery charging regulator IC1 (pin 5) through Q1. Biasing for Q1 is accomplished by D2, R1 and R2 when AC power is applied. When running on battery power Q1 is biased off by R1, R2 and D3, this prevents the battery from trying to power the battery charger regulator and IC2 that informs the monitor of the loss of AC.

The output of switching regulator IC1 pin4 is rectified and filtered by D4, C4 and L1 then fed to the battery through current sense resistor R3 and fuse F302 to J302 pin 1 (VBAT+). The battery float charge voltage is maintained at 13.2 volts except for fast charge which is regulated at 14.4 volts. The output is also monitored for over current conditions. These parameters are controlled by IC3 and associated circuitry. When the battery charge current exceeds 120mA of current IC3 pin 7 goes high which biases Q2 on, this in turn shorts out R12 which affects the feedback control (FB) to IC1 (pin 1). With R12 shorted out the control resistors R14 and R13 set the output voltage to 14.4 volts. When the charge current lowers IC3 pin 7 goes low which biases Q2 off, this puts R12 back into the feedback control which now consists of R12, R13 and R14 setting the voltage to 13.2 volts. When more than 600mA of current flows through R3, IC3 pin 1 shorts IC1 pin 2 to ground which shuts IC1 off until its next switching cycle, when the current reaches a safe level IC3 pin 1 allows IC1 to remain on.

The voltage switched by Q1 is also fed to IC2 as VCH (Voltage Charge). The output of this 5 volt regulator provides the LINEST (Line Status) signal to the main board at E302 pin 3. With AC applied, LINEST is high. LINEST goes low when the AC is disconnected. The LINEST line is also routed to the power on/off circuitry. See *Power On/Off Control Circuitry* on page 19.

Battery Operation

9.1.3

Without AC power there will be no secondary voltage rectified by D1. Power for the monitor will be supplied from the battery at J302 pins 1 (VBAT+) and 2

(VBAT-). The battery power will conduct through D3 and F301 to VIN at E302 pin 1 to the 2471 main board. The arrangement of R1, R2 and D2 bias Q1 off in this condition which prevents power from reaching IC1, IC2 and IC3. The output of IC1 is also protected by D5 which is now reverse biased, the bridge D1 is also reverse biased and prevents T301 from discharging the battery. With no voltage at IC2 the LINEST will be low which indicates to the main board that there is no AC power.

2726 Power Supply Board

9.2

The 2726 Power Supply Board contains the circuitry needed to power the monitor from the AC Mains (line voltage). The power supply board also connects to the battery and contains the battery charging circuitry.

AC Mains and Battery Operation Overview

9.2.1

The Model 520A can be powered from its internal 12 volt battery or from the AC Mains. The green (sine wave icon) front panel indicator illuminates when the line cord is connected and the rear panel power switch is in the “I” (On) position. This indicates that AC Mains power is reaching the power supply, that the battery is being charged, and that if the monitor is turned on, it is being powered from the line voltage.

If AC Mains power is removed by unplugging the line cord or setting the rear panel power switch to the “0” (Off) position, the monitor will operate for up to four hours from its internal 12 volt lead-acid battery. As the battery voltage runs low (<11.5 volts), the red battery indicator on the front panel illuminates. At this point, the AC Mains should be reconnected to power the monitor and charge the battery.

If the monitor continues to be powered from a battery in a low voltage situation, at approximately 11 volts, a continuous alarm sounds for thirty seconds while the Message Center displays BATTERY EXHAUSTED CONNECT LINE CORD. If this alarm/message is ignored, the monitor displays will shut down and the battery indicator will flash on and off about every 5 seconds. If AC power is now restored, the monitor will re-initialize (go through the power up and self-test routines) and resume normal operation. However, continued battery operation will eventually activate a hardware low voltage circuit (<8.5 volts) that shuts the monitor off to prevent damage to the battery. Once the unit is shut down with the hardware circuit, the AC Mains must be connected and the front panel POWER key pressed before the monitor will turn back on.

AC Mains Operation

9.2.2

The AC Mains voltage enters the monitor at the rear panel Power Entry Module (PEM). This device contains a built in RFI power line filter, a double-pole single-throw switch that opens and closes both AC input lines, fuses, and an input voltage selection card.

The filtered, switched and fused output of the Power Entry Module is fed to the primary coils of the rear panel mounted system transformer, T301. The secondary output from T301 is rectified by D1 (bridge rectifier) and filtered by C1. The loaded DC voltage at this point is approximately 20 volts.

The 20 DC volts is fed to the main board through fuse F301 to connector E302, and is switched to the battery charging regulator IC1 (pin 5) through Q1B. Biasing for Q1B is accomplished by C2, R1, R2, C8 and Q1A when AC power is applied. When running on battery power Q1B is biased off by R2 and Q1A, this prevents the battery from trying to power the battery charger regulator and IC2 that informs the monitor of the loss of AC.

The output of switching regulator IC1 pin4 is rectified and filtered by D4, C4 and L1 then fed to the battery through current sense resistor R3 and fuse F302 to VBAT+ (J302 pin 1). The battery float charge voltage is maintained at 13.2 volts except for fast charge which is regulated at 14.4 volts. The output is also monitored for over current conditions. These parameters are controlled by IC3 and associated circuitry. When the battery charge current exceeds 120mA of current IC3 pin 7 goes high which biases Q2 on, this in turn shorts out R12 which affects the feedback control (FB) to IC1 (pin 1). With R12 shorted out the control resistors R14 and R13 set the output voltage to 14.4 volts. When the charge current lowers IC3 pin 7 goes low which biases Q2 off, this puts R12 back into the feedback control which now consists of R12, R13 and R14 setting the voltage to 13.2 volts. When more than 600mA of current flows through R3, IC3 pin 1 shorts IC1 pin 2 to ground which shuts IC1 off until its next switching cycle, when the current reaches a safe level IC3 pin 1 allows IC1 to remain on.

The voltage rectified by D1, D2 and filtered by C2 is fed to IC2 as VCH (Voltage Charge). The output of this 5 volt regulator provides the LINEST (Line Status) signal to the main board at E302 pin 3. With AC applied, LINEST is high. LINEST goes low when the AC is disconnected. The LINEST line is also routed to the power on/off circuitry.

Battery Operation

9.2.3

Without AC power there will be no secondary voltage rectified by D1. Power for the monitor will be supplied from the battery at J302 pins 1 (VBAT+) and 2 (VBAT-). The battery power will conduct through D3 and F301 to VIN at E302 pin 1 to the main board. R2 and Q1A bias Q1B off in this condition which prevents power from reaching IC1, IC2 and IC3. The output of IC1 is also protected by D5 which is now reverse biased, the bridge D1 is also reverse biased and prevents T301 from discharging the battery. With no voltage at IC2 the LINEST will be low which indicates to the main board that there is no AC power.

2471 & 2775 Main Board

9.3

The Main Board contains all the analog and digital circuitry that controls the sensor, external communication and front panel display. The isolated power supplies, microprocessor circuits and memory are all contained on this board. The earlier 2471 Main Board and newer 2775 Main Board are identical with exceptions noted in the text.

Power On/Off Control Circuitry

9.3.1

The Model 520A power on/off control circuitry consists of the VBACK supply (regulated by IC12), IC10, IC11 and the **POWER** key. (See page 4 2471 schematic, page 4 on 2775 schematic.)

When the battery or AC Mains is first applied to the power supply board (via VIN J102 pin 1), VBACK goes to +5 volts, provides power to IC10 and IC11, and through the C26 and RP4 (pins 3,4) network at IC10 pin 8, sets IC10 pin 2 to a logic Low.

The ON/OFF line is brought Low each time the front panel **POWER** key is pressed. This sends the output at IC11 pin10 High. This Low-to-High transition clocks the (#1) D flip-flop portion of IC10. The /Q1 output at IC10 pin 2 goes High and with each successive press of the **POWER** key, this output toggles to the opposite level (Low or High). A High turns the Model 520A on and a Low shuts it off.

While the /Q1 output at IC10 pin 2 is High, the MOSFET Q7 is turned on and pulls the gate of MOSFET Q8 to ground, thus causing Q8 to conduct as well. With Q8 conducting, the currently active monitor power source—either the AC Mains derived supply or the battery supply will flow through Q8 to the voltage input (pin 7) of the Pulse Width Modulator IC9. The output IC9 pin 6 will oscillate (at the frequency set by R13 and C15). This causes Q5 to switch on and off and provide a path to ground through the primary coils of T1 for the supply (Mains or battery) at T1 pin 12. Current flowing in the primary is measured at IC9 pin 3 and the duty cycle of the pin 6 output will vary with the load on the transformer.

Current flow in the transformer primary induces current in the three secondary coils and creates the ± 12 volt analog supplies (+V12 and -V12), the VRAWI that powers the isolated RS232 circuitry, and the +5 volt VCC supply that power the remaining circuits in the monitor. The Model 520A turns on. The +V12 and -V12 supplies are rectified and filtered by D2, D4, C10, C11, C12 and C13. The +V12 is regulated by IC7 and the -V12 by IC8. The VDD supply is rectified by D3, filtered by L1, C9 and C20 and fused at F1, and in addition, a feedback loop to IC9 contains VR1 which is factory adjusted to produce a +5.00 volt ± 0.05 volt VCC supply (measured under load).

Once the monitor powers up, a SYNC signal toggles Q9 on and off causing a timing pulse to be transmitted across C19 and C15 to the input at IC9 pin 4. This has the

effect of synchronizing the output of the pulse width modulator with the data sampling operations of the analog board and keeps power supply switching spikes from interfering with those operations.

Power Supplies

9.3.2

The secondary pins 7, 8, 9 of T1 form a center tap transformer, the voltage is rectified by D2 and D4 then filtered by C10, C11, C12 and C13. The dual 12 volt supplies, +V12 and -V12 which are generated from this voltage are regulated by IC7 and IC8 respectively. The secondary winding of pins 5 and 6 of T1 are rectified by D3 and filtered by C9, L1 and C20, this voltage designated as VCC (+5 volts) acts as reference for IC9, supplies power for the opto isolator non-isolated side and powers other circuitry on the board.

The secondary winding consisting of pins 2 and 3 are rectified and filtered by D1 and C1. The rectified voltage at this point is approximately 7 volts DC and is regulated to 5 volts by IC2. This isolated supply powers the isolated portion of the opto-isolators and the RS232 driver chip IC1. The unregulated voltage VRAWI is sent to the rear panel connector J101.

The backup voltage (VBACK) is regulated by IC12 from the VIN supply. Capacitors C22 and C27 serve as filters and D17 allows VCC to power VBACK circuitry when the monitor is on. At this point D18 is biased off so IC12 is idle. When the monitor is turned off and VCC collapses D18 is then forward biased and IC12 now supplies VBACK circuitry, D17 at this time is reverse biased and prevents power from reaching VCC.

The saturation sensor LEDs derive their power (LEDSRC) from the current regulator IC32. (See sheet 3 on 2471 schematic, page 2 on 2775 schematic.) Resistor R31 limits the maximum current draw to 45 mA (nominal draw 35 mA). Regulator output is filtered by C85 and L2. The charge stored on capacitor C1 supplies the 290-350 mA peak currents that can occur when the sensor LEDs are turned on. Diode D12 prevents the regulator output from exceeding +7.5 volts while the fuse, F2, provides current limit protection in the event of a regulator circuitry failure. The RP10 (pins 1,8 and 5,6) divider network provides the CPU (via IC33) with a means to monitor the LEDSRC status.

The +VA and -V12 supplies are regulated to +V5 (+5 volts) and -V5 (-5 volts) by IC44 and IC43 respectively. These supplies are used by the 20 bit ADCs, the 8 bit ADC and other circuits associated with them.

Voltage References

9.3.3

A +2.5 volt precision reference supply, generated by IC35 from the 12 volt +V12 supply, is used as a reference voltage for the ADC chips IC37 and IC34.

The +2.5 volt output from IC35 pin 6 is fed to the non-inverting input of amplifier IC36 pin 3. Resistors R32 and R33 combine for a gain of 1.617 that provides a +4.096 volt reference (approximately) supply, 4VREF, at IC31 pin 1.

The 4VREF is fed to IC30 pin 13 which is set up as a unity gain inverting buffer amplifier, therefore the output at IC30 pin 14 is -4 volts. This -4 volts is used by IC29 as a reference voltage for VLED (Voltage LED) and CNTRST (Contrast) controls. (See sheet 4 on 2471 schematic, sheet 3 on 2775 schematic.)

Preserving RAM and Real Time Clock Data

9.3.4

The NAND gate output at IC11 pin11 will be Low when the monitor is on (IC10 pin 2 is High) and High when the monitor is off (IC10 pin 2 is Low). (See sheet 4 of 5 on schematic.) This PWRON* (Power On) signal is used to prevent corruption of RAM and clock data when the monitor is turned off. It does this by going High and therefore denying CPU access to the RAM and clock chips so that as the power supplies crash when the monitor is turned off, the CPU cannot send erroneous data to the chips.

Whenever the CPU is writing information to the RAM or Real Time Clock chips, the CPU momentarily sends the OFFDIS (Off Disable) line High. The High going level appears at IC11 pin 2. Since the monitor is powered on, IC11 pin1 will be High. This means IC11 pin3 momentarily goes Low, Q6 starts to conduct and IC10 pin 4 goes High. In this reset condition the /Q1 output #1 flip-flop of IC10 (pin 2) will be held high even if the user presses the **POWER** key and clocks the flip-flop. In effect, the CPU is not allowing the monitor to be turned off. The Low at IC11 pin3 will last for the duration of the RC time constant set by C21 and RP4 (pins 5,6). These values were chosen to produce a time-out longer than the time necessary to complete the write to RAM or Real Time Clock operation. After the RC time-out, IC10 pin 4 returns Low and a press of the **POWER** key will toggle the /Q1 output of IC10 and the monitor will turn off. This is done to prevent corruption of RAM and Real Time Clock data.

Low Battery Voltage Shutdown

9.3.5

The CPU monitors the battery voltage and provides the user with a low battery indicator (E3), messages and alarms. However, if these are ignored, a hardware circuit will take over and shut off the monitor before the battery is damaged.

The pulse width modulator IC9 requires at least 7.6 volts at pin 7, its voltage supply, in order to operate. This pin typically draws 10 mA of current. The resistance of the R63 and Q8 combination is approximately 114 ohms. This equates to a voltage drop of approximately 1.14 volts. Therefore if the battery voltage drops under 9.0 volts (approximately), IC9 will not have sufficient voltage to operate and will shut down. Shutdown of IC9 stops current flow through transformer T1 and the secondary supplies shut down, effectively turning off the monitor.

When IC9 shuts down, its VREF output at pin 8 is pulled Low. This forward biases D6 and causes the NAND gate output at IC11 pin4 to go High. The #2 flip-flop of IC10 is clocked, and the High at the D2 input (because /Q1 is High) is transferred to the Q2 output at pin 13. The High at pin 13 Sets the #1 flip-flop causing the /Q1 pin 2 output to go Low. This Low shuts off both MOSFETs Q7 and Q8, thereby blocking any supply voltage from IC9 pin 7. Normally, pressing the front panel **POWER** key would clock flip-flop #1 (at pin 3) and return the pin 2 output High—but the High output at pin 13 keeps the #1 flip-flop Set—and the **POWER** key has no effect.

If at this point the AC MAINS is reconnected, MOSFET Q8 continues to block current from IC9 pin 7 and the monitor remains off. Connecting the AC Mains does however send the LINEST signal High. This High Line Status signal is brought to IC10 pin10 where it Resets the #2 flip-flop, sending IC10 pin 13 Low and removing the Set condition from flip-flop #1. Now, if the front panel **POWER** key is pressed, flip-flop #1 is clocked, IC11 pin11 goes High, MOSFETs Q7 and Q8 turn on, the supply to IC9 pin 7 is restored, the pulse width modulator restarts, energizes T1, and the monitor turns back on.

Timing Sequencer

9.3.6

A 14 stage divider, IC39, acts as a timing sequencer. (See sheet 3 on 2471 schematic, sheet 2 on 2775 schematic.) A 3.276 MHz crystal Y2, provides a Clock Sequence (CLKSEQ) to IC39 pin 10. The IC39 pin 11 RESET input line resets IC39 on monitor power up.

The IC39 pin 9 Q1,output provides a clock input signal to the tone generator board located at IC27. The IC39 pin 3 Q14 output provides a 5 ms interrupt (INT5MS) for IC18. The Q4-Q11 outputs of IC39 become inputs to the Data Sampling Controller IC42.

Data Sampling Controller

9.3.7

The IC39 Timing Sequencer's Q4-Q11 outputs become inputs to IC42, a PEEL (Programmable Electrically Erasable Logic) device. The PEEL uses the CLK and D0-D6 inputs, and the SC1 and SC2 inputs, to control data sampling by providing sensor LED drive signals and demultiplexing for the signals returning from the saturation sensor.

The waveforms in Figure 3.(with the exception of CLK) are only valid when both the SC1 and SC2 inputs are low. The System Calibration inputs (SC1 and SC2) generated by the microprocessor, are kept low, except that they are toggled high/

low, during a Probe Off Patient alert, and during a system power up self-test. See *Calibrating the 20-Bit Analog-to-Digital Convertors* on page 25.




Figure 3. Figure 1. Front-End Timing

The data sampling sequence consists of:

- a. turn on the Red LED (RDLED*) and the Analog Sample line (ASAMP*)
- c. sample the Red LED return signal (RDSMP*)
- e. turn on the Infrared LED (IRLED*) and the Analog Sample line (ASAMP*)
- f. sample the Infrared LED return signal (ISMP*)
- h. turn off the Infrared LED, the Analog Sample line, and stop sampling
- i. repeat the process starting at step a.

The Analog Sample (ASAMP*) line is used to nullify the effects of any ambient light signals returning from the sensor. See *Sensor Photodiode Return Path* on page 24.

The IC42 INSIG* and SIGND* outputs are used in conjunction with the SC1 and SC2 inputs. See *Calibrating the 20-Bit Analog-to-Digital Convertors* on page 25.

The IC42 pin 14 external sequencer (SYNC) line is equivalent to the PEEL's D1 input. It provides a "sync" pulse to the pulse width modulator on the power supply board in order to keep power supply switching spikes from interfering with data sampling operations.

Sensor LED Drive Circuits

9.3.8

The VLED line voltage is derived from IC30 pin 8 which is controlled by the Digital to Analog Converter IC29. (See sheet 4 on 2471 schematic, sheet 3 on 2775 schematic.) When address line A0 is high (IC29 pin 6) and both WR* (IC29 pin 16) and DACCS* (IC29 pin 15) are low the D/A Converter is enabled. The data on lines

D0-D7 (IC30 pins 14-7) now control the output voltage of IC30 pin 8 (VLED) based on the VREFB voltage on IC29 pin 18.

When the RDLED* signal at IC42 pin 18 goes low (logic 0), Q14 turns off and the Red LED signal (VLED from IC30 pin 8) at R37 is divided by R37 and R41, finally causing IC36 pin8 to go high. This positive output turns Q11 on and current flows from the LED source (LEDSRC), through the Red LED in the sensor (it turns on) returning as LED1SK (LED1 sink) across Q11 and the current limiting resistor R53 to ground. (See sheet 3 on 2471 schematic, sheet 2 on 2775 schematic.)

When RDLED* returns high (logic 1), Q14 is biased on, forcing IC36 pin8 to ground potential: Q11 is biased off, and as a result, the Red LED in the sensor is also off.

The Infrared LED drive circuit operates in the same manner as the Red LED drive discussed above. The IRLED* signal at IC42 pin 17 activates Q13 the LED2S signal causes a positive signal at IC36 pin 14, and current can flow from the LEDSRC supply through the sensor's Infrared LED, Q12, and the current limiting R52 to ground.

Sensor Photodiode Return Path

9.3.9

Light, from the sensor's Red or Infrared LED, shines through the pulsating vascular bed (the patient's finger, toe, etc.) placed between the LEDs and the photodiode. Some of this light emerges from the tissue and impinges on the photodiode, causing the photodiode to conduct current. IC40 pins 1-3 is set up as a differential amplifier that converts this input current to a voltage at the amplifier output. The sensors are wired such that photodiode current produces a positive voltage at IC40 pin 1¹.

The voltage at IC40 pin 1 is presented to an analog switch IC41 pin 6. This switch is controlled at pin8 by INSIG* (Input Signal) from IC42, and will be closed (IC41 pins6 and7 connected) except if the monitor is in a Probe Off Patient condition or is undergoing its Self-Test at system power up. The switch IC41 pins 9-11, controlled from SIGND* (Signal Ground) at IC42 will be open (no connection between IC41 pins10 and 11) except as noted above for the switch at IC41 pins 6-8. As a result, the IC40 pin 1 voltage passes undisturbed to the high pass filter consisting of R59 and C90.

As shown in Figure 3., the ASAMP* signal is active whenever either sensor LED is turned on. This causes Q15 to turn off and the charge at C90 passes through the unity gain buffer to IC40 pin 7.

If the signal at IC40 pin 7 is the product of the Red LED being turned on, then RDSMP* from IC42 pin 12 will go low and close the switch at IC41 pins 2-3, thereby presenting the signal to a sample and hold circuit consisting of R54 and

1. The Model 520A uses SuperBright™ sensors. If a non-SuperBright™ sensor is connected, IC40 pin 1 will go negative—a condition that causes an “Incompatible Probe” message to be displayed.

C100 (that maintains the signal until next sample pulse arrives), a gain stage, (IC38 pin 1), a filter/divider network (C87, R45 and R46), and finally, to the Red channel Analog-to-Digital Convertor (ADC) IC34.

If the signal at IC40 pin 7 is the product of the Infrared LED being turned on, then ISMP* from IC42 pin 13 will go low and close the switch at IC41 pins 14-15, thereby presenting the signal to a sample and hold circuit consisting of R55 and C96 (that maintains the signal until next sample pulse arrives), a gain stage, (IC38 pin 7), a filter/divider network (C88, R49 and R50), and finally, to the Infrared channel Analog-to-Digital Convertor IC37.

Again referencing Figure 3., the ASAMP* line returns to a logic high when neither LED is being driven, causing Q15 to turn on. With Q15 conducting, any charge at C90 is discharged to ground and the next pulse will charge C90 from a known level. If it were not for Q15, any charge remaining on C90 from the previous pulse or from ambient light reaching the photodiode would be added to the charge from a new pulse—creating measurement errors.

Calibrating the 20-Bit Analog-to-Digital Convertors

9.3.10

The 20-bit ADCs are calibrated as part of the system self-test which occurs each time the monitor is turned on. At power up, the microprocessor sets the CAL line high. The System Calibrations input SC1 is set high and SC2 is reset to a logic low. The CS5503 ADC will not operate while the CAL line is high. On the falling edge of the CAL signal, the ADC will initiate a calibration cycle determined by the state of the SC1 and SC2 inputs.

The high at SC1 and the low at SC2 cause the Data Sampling Controller, IC42, to set INSIG* high and reset SIGND* to a logic low. The high INSIG* opens the switch at IC41 pin8 so that IC41 pins 6 and 7 are no longer connected—disconnecting the returning photodiode signal from the rest of the circuitry. The low SIGND* signal closes the switch at IC41 pin9 and as a result, the input to the C90-R59 high pass filter (and thus the entire ADC input circuitry) is brought to ground potential.

The CAL line (which went high at power up) is reset low and ADCs IC34 and IC37 begin their calibration cycles. Because the analog input circuitry is grounded via SIGND*, only circuit offset voltages can be present at the (pin 9 AIN) input. The calibration cycle sets the ADC “zero” point to equal this voltage, thus compensating for any circuitry offsets. The ADC then sets its “full scale” point to equal the voltage at its VREF (pin 10) input. This completes the calibration cycle.

The ADC can now start sampling its input and converting it to a 20-bit digital word. The processor resets SC1 to a logic low, causing IC41 pin9 to open and IC41 pin8 to close. The photodiode signal can now reach the ADCs. See *Sensor Photodiode Return Path* on page 24.

20-Bit Analog-to-Digital Conversion

9.3.11

Data from the Red and Infrared channels is sampled by the 20-bit measurement ADCs, IC34 and IC37 respectively. The analog input at pin 9 is converted to a digital representation with 20-bit resolution based on the input magnitude.

The CS5503 convertor continuously samples its input, converts the value to a digital word, puts the word in its output buffer (overwriting previous buffer contents), then repeats the process by again sampling its input. The frequency of the sample/convert/overwrite-buffer sequence is based on the 3.2768 MHz clock signal at the ADC pin 3 (CLKSEQ) input.

The microprocessor starts a read cycle of the Infrared channel by bringing IC37 pin 16 (Chip Select Channel 1) low. A Red channel read starts when IC34 pin 16 (Chip Select Channel 2) is brought low.

On the falling edge of the ADC's CS*, the output word's MSB (most significant bit) appears at the pin-20 SDATA (Serial Data) output. The SDATA line connects directly to the microprocessor's serial input (RXS) pin. The remaining bits (in descending order) are output from SDATA with subsequent falling edges of the Serial Clock (SCLK) input at pin 19. The SDATA output automatically goes to a 3-state (high impedance) condition after completing a word transmission, thus freeing the data line for other uses (i.e., the other ADC channel).

The Serial Clock speed is controlled through the digital board PEEL IC18. This clock rate is significantly slower than the ADC sampling rate. As a result, the ADC rewrites its output buffer with new information at a faster rate than the data can be read from the buffer. No conflict occurs, however, because while CS* is low (during the read cycle), the ADC does not update its output buffer—the current word is not overwritten. After the processor receives the entire word, it allows the convertor's CS* to return high, and the ADC resumes its sample/convert/overwrite-buffer cycle.

Sensor Status Decoding and Conversion

9.3.12

The microprocessor monitors several sensor parameters in addition to the Red and Infrared data channels. It monitors the status parameters, as well as the voltage of the monitor's internal battery.

The 8-to-1 multiplexor, IC33, decodes the A0AUX-A2AUX input address lines and connects one of eight status parameter inputs (labeled channels 0-7 at IC33) to the multiplexor output at IC33 pin 3. Resistor R29 and diode D13 prevent negative voltages from reaching the input to the analog-to-digital convertor, IC31.

IC31 is an 8-bit analog-to-digital convertor with a serial data output. While the IC31 Chip Select (\overline{CS}) input is high, the CLK input and DOUT output are in 3-state mode. When \overline{CS} is brought low (under processor control), the most significant bit (D7) of the *PREVIOUS* data conversion becomes available at the DOUT pin. The

remaining bits (D6-D0) are shifted out on subsequent falling edges of the CLK input. On the clock pulse following the one that shifts out the least significant bit (D0), the CLK and DOUT lines are returned to 3-state and the ADC performs a new conversion based on the input it receives from the IC33 channel selected by the A0AUX-A2AUX input address lines.

The ADC sample/convert/store-result cycle is based on internal chip timing and not the CLK input which (along with \overline{CS}) only controls serial data output. Thus the \overline{CS} line is free to return high once the ADC cycle begins.

Sensor Status Parameters

9.3.13

The sensor (and battery) status parameters input to the multiplexor IC33 are described below. Note that channel number refers not to the IC33 pin number, but to the signal label (e.g., channel I0 signal resides at IC33 pin 13).

Channel I0: ADCV12.

This is an extra input to the multiplexer IC33 pin 13. It is unused as of this writing.

Channel I1: Auxiliary Input.

This is an extra input to the multiplexer IC33 pin 14. It is unused as of this writing.

Channel I2: Photodiode DC Level.

Resistors R40, RP8 (pins 1, 2) and capacitor C97 form a voltage divider and low pass filter that provide a measure of the mean DC level at the output of the photodiode current-to-voltage amplifier IC40 pin 1. This channel (IC33 pin 15) is used in determining ambient light interference. If this line is examined while the sensor's Red and Infrared LEDs are turned off, then any DC level at IC40 pin 1 must be the result of ambient light impinging on the photodiode. If the DC shift is in excess of limits set in the software, a Light Interference message appears on the monitor's display.

Channel I3: Sensor LED Supply Voltage.

This channel, at IC33 pin 12, monitors the sensor LED supply voltage through a voltage divider consisting of RP10 (pins 5, 6) and RP8 (pins 1, 8). If a fault occurs that causes the LED supply fuse F2 to blow, or if the sensor wires are shorted, this channel reports the condition and the Message Center displays PROBE FAILURE.

Channel I4: Incompatible Probe Detection

The input at IC33 pin 1 provides the processor with an incompatible probe indicator (words probe and sensor are interchangeable). The photodiode output voltage at IC40 pin 1 will be positive if a SuperBright™ sensor is connected to the monitor. This positive signal passes through the high pass filter of C73 and RP10 (pins 7, 8) to the amplifier inverting input IC36 pin 6, where it is summed with the -5 volt output from IC43. The resultant voltage at IC33 pin 1 will be approximately +2.81 volts with no photodiode input. This voltage drops as the (positive) photodiode signal at IC40 pin 1 increases. If the IC40 pin 1 signal goes negative, as would happen if a non-SuperBright™ 8600 series sensor were connected to the monitor,

the cumulative effect of the -5 volts and the negative photodiode signal passing through D14 and R34 would be to send the IC36 pin7 output to its positive rail (+12 volts). The processor recognizes this higher voltage and causes an INCOMPATIBLE PROBE message to be displayed in the Message Center.

Channel I5: Infrared LED Cathode Voltage.

A low pass filter consisting of RP10 (pins 3, 4), RP8 (pins 1, 7) and C89 provides a means to measure the cathode voltage of the sensor's Infrared LED. If the channel at IC33 pin 5 is sampled the monitor can determine if the LED is open circuit (zero volts at IC33 pin 5) or operational (approximately 2.5 volts at IC33 pin 5). If not operational, display messages such as PROBE I.R. LED FAILED or ERROR - FAULTY PROBE are generated.

Channel I6: Battery Supply Voltage

The monitor's internal battery voltage is divided down by RP3 (pins 7, 8) and RP8 (pins 1, 5). The voltage at IC33 pin 2 is monitored and if its magnitude is less than a predetermined value (encoded in the software) the monitor lights and/or flashes its front panel battery indicator. This provides the user with a low battery warning.

Channel I7: Red LED Cathode Voltage.

A low pass filter consisting of RP10 (pins 1, 2), RP8 (pins 1, 8), and C92 provides a means to measure the cathode voltage of the sensor's Red LED. If the channel at IC33 pin 4 is sampled the monitor can determine if the LED is open circuit (zero volts at IC33 pin 4) or operational (approximately 2.5 volts at IC33 pin 4). If not operational, display messages such as PROBE RED LED FAILED or ERROR - FAULTY PROBE are generated.

Microprocessor

9.3.14

A Hitachi HD64180RP microprocessor directs the actions of the Model 520A Pulse Oximeter. (See sheet 2 on 2471 schematic, sheet 1 on 2775 schematic.) The processor, IC16, is operated at 6.144 MHz (half the 12.288 MHz frequency of crystal Y1), has an 8-bit data bus and a 19-bit address bus (the 520A uses only 18-bits). The microprocessor also provides two asynchronous serial communication channels, a clocked serial I/O port and various interrupt and control signals. The +5 volt VCC supply to the processor is first sent through inductor FB1, a ferrite bead, before powering the chip at IC16 pin 32.

Memory

9.3.15

The Model 520A system software is located at IC17, a 27C010 EPROM. The 32 K byte RAM, IC20, stores trend data, system power up settings (averaging times, serial output parameters, etc.), and provides an area for system (scratchpad) memory requirements. Since IC20 is powered from the VBACK supply, RAM memory is retained when the monitor is turned off and it becomes available again when the monitor is turned back on.

The ROM at IC17 is read when its Chip Enable line (IC17 pin 22) is brought low by the ROMCS* signal at IC25 pin 3, and the processor brings its Read line (IC16 pin 63) low—thereby activating the ROM Output Enable line at IC17 pin 24. Under these conditions, ROM data from the specified address bus location is made available to the data bus for use by the processor.

The RAM (IC20) is activated when its Chip Select line RAMCS* (IC20 pin 20) is brought low, via IC25 pin 8. When the ME* line at IC25 pin 5 is low, and the inverse of address line A17 at IC25 pin 4 is low, output pin 6 of IC25 will go low. This in turn will drive IC25 pin 10 low, with PWRON* low at IC25 pin 9, IC25 pin 8 will be low (RAMCS*). If at that time, \overline{OE} (IC20 pin 22) is low, a RAM Read occurs, whereas a RAM Write will occur if \overline{WE} (IC20 pin 27) is low.

Real Time Clock (RTC)

9.3.16

A Real Time Clock provides the Model 520A the ability to time stamp collected (printed) trend data. (See sheet 4 on 2471 schematic, sheet 3 on 2775 schematic.) The 32.768 kHz crystal, Y3, provides the timing signals for the clock chip, IC24, which is powered from the VBACK supply so that the clock can continue to keep time when the monitor is turned off (provided the monitor's 12 volt internal battery is connected and maintains at least a nominal charge).

The RTC is activated when its Chip Select line (CS0*) at IC24 pin 2 is brought low. With the monitor on PWRON* will be low (IC25 pin 13), the RTC* line will be brought low by the processor through IC22 pin 14, these lines drive RTCCS (Real Time Clock Chip Select IC25 pin 11) low. If at this time, the RD* (IC24 pin 8) input is low, a RTC Read occurs, whereas a RTC Write will occur if the WR* (IC24 pin 10) input is low. Addressing is handled by A0-A3 (pins 4-7) and data I/O through D0-D3 (pins 14-11).

Sound generator

9.3.17

The programmable tone controller IC27 is replaced with the Tone Generator Replacement board 2755-01, this is used to drive the monitor's audio circuit. The tone generator is clocked by IC39 pin9 from the 1.638 MHz signal, Q1. The tone generator is enabled by the processor when IC22 pin12 is brought low. While CE* is low, WR* is brought low and data bus information including frequency (pitch) and attenuation (volume) is accepted by the tone generator. The Ready signal (IC27 pin 4) goes low while accepting data and the processor is put into a Wait state until IC27 finishes its task; then Ready returns high and the processor continues its operations.

The AUDOUT output at IC27 pin 7 drives the audio amplifier IC26. The amplifier output is coupled through capacitors C55 and brought to J109 as the SNDOUT (Sound Out) line. The speaker which is mounted in the chassis is connected to J109, LS1 is not installed on the 2471 Main Board.

Keypanel Interface**9.3.18**

The 12 keys (switches) on the Model 520A front panel are connected to the 2471 Main Board through a ribbon cable at J104. Each key (except **POWER**) is connected to an 8-bit latch (either IC14 or IC15). When any of these keys is pressed, the corresponding latch input is brought low. The processor continually reads the status of these latch outputs, the RDKEY* enables IC14 when low and the RDKEY_2* line enables IC15 when brought low.

The **POWER** key ON/OFF signal is sent through J104 pin 15. The AC Line Status signal, LINEST, is generated by the power supply, and is high (+5 volts) when the monitor is connected to the AC Line (Mains) and the rear panel switch is set to "I". This +5 volt level is sent to LED D2 (the green ~ indicator) on the keypanel via J104 pin 16. The LINEST signal is also input to the latch at IC14 pin 8 so that the CPU can detect if the unit is operating on AC line power (IC14 pin 8 high) or on battery power (IC14 pin 8 low). Diode D9 keeps current from backfeeding into IC14 when the monitor is turned off but still connected to the AC line.


The 2MIN (D3), AUDIO OFF (D4), ALERT (D5), and LOW BAT (D1) LEDs are driven by the 8-bit latch IC13. When each of the corresponding signals is driven high the appropriate LED on the keypanel is illuminated. See *I/O Device Controller* on page 31.

Display Interface**9.3.19**

The display is connected to the 2471 Main Board at J107. It is controlled by the processor using the RD* (Read), WR* (Write), and DISPCS* (Display Chip Select) lines. Data bits D0-D7 are used as input/output lines and A0 is used in conjunction with the RD* and WR* lines to distinguish between read and write operations as listed below.

A0	RD* Low	WR* Low
High	Status Read	Command Write
Low	Data Read	Data Write

Table 1. Display Command/Data table.

The CNTRST (Contrast) line is controlled by pressing the front panel  key. When depressed and held the processor controls the digital to analog converter IC29, line A0 is brought low along with WR* and DACCS*, the data on D0-D7 controls the voltage at IC30 pin 7 which can vary from 0-4 volts. The output of IC30 pin 7, along with 4VREF feeds a summing amplifier (IC30 pins 1, 2, 3). The output of the summing amplifier IC30 pin 3 controls the base of Q10 which in turn controls the contrast of the display through a variable negative voltage.

The backlight for the display is controlled by the DSPBR (Display Bright) line. When DSPBR is high the gate of Q16 is biased off, current flows from Vcc through R51 to IC45. This sets the backlight for low illumination. The illumination of the backlight is made greater when DSPBR is made low, this biases Q16 on which essentially shorts out R51 allowing more current to flow into IC45 creating a full backlight.

I/O Device Controller


9.3.20


The A/D Converter Chip Selects, serial A/D Chip Selects, Sensor Status Decoding and NEXT* line are all controlled by IC28 when selected by the OPORT line (IC21 pin 10). The OPORT line will go high when the L1* and WR* line both go low at IC23 pins 13 and 12, this will send output pin IC23 pin 11 low which drives inverter IC21 pin 10 high enabling IC28.


A 3 to 8 decoder is used to control the DACCS*, RTC*, DISPC*, AUD*, KEYS*, L1*, L2*, 2KEYS* lines. when the IOE* line goes low and the LIR* line goes high being inverted by IC21 pin 2 and presented to IC22 pin 5 as a low enable line IC22 is enabled, Q0-Q7 will be driven low depending upon the A4, A5 and A6 lines on pins 1, 2, 3 respectively on IC22.

With the LPORT line high IC13 is enabled, this latches the data on lines D0-D7 (1D-8D pins 2-9) on its output pins 19-12 (1Q-8Q respectively), the outputs correspond to the following eight lines:

CAL-used by the A/D Convertors on power up to compensate for front end voltage offsets.

TML-Two Minute LED drives the  LED on the front display.

AOL-Audio Out LED drives the  LED on the front panel.

ARL-Alert LED drives the  LED on the front panel.

KJL*-drives Q17 when high which in turn drives the Alert Bar LEDs via J105.

BTL-Battery Low  LED on the front panel.

DSPBR-(not used on the 520A)

OFFDIS-sent to the Power On/Off section of the circuitry to prevent the monitor from being turned off while writing to RAM.

Watchdog Timer

9.3.21

The Watchdog Timer provides a system reset function in the event a hardware or software “glitch” occurs. The PEEL IC18 forms the heart of the Watchdog circuit. (See sheet 2 on 2471 schematic, sheet 1 on 2775 schematic.)

At powerup and at specific intervals thereafter, the microprocessor outputs a logic high to IC18 pin 8, WDOG (Watchdog). The WDOG signal combines with other signals within the PEEL and as a result the Watchdog Clear (WDCLR) open drain output at IC18 pin 12 is continually brought low. This discharges the capacitor C46 before it can charge up (via RP2 pins 1 and 2) past the input threshold of IC21 pin 9.

If the WDOG pulse does not appear at regular intervals, as the result of a software or hardware problem, the RC charges up and IC21 pin 8 goes low producing a Watchdog Time-out (WDTO*) input at IC18 pin 11. WDTO* combines with other signals within the PEEL and causes the open drain Master Reset (MR) output at IC18 pin 13 to be brought low.

A low MR causes C45 to discharge, forcing IC21 pin6 high. This sends a reset pulse to the system. It also sets the Reset Input (RESIN) signal at IC18 pin 9 high which causes RESET* at IC18 pin 19 to activate low. The active RESET line causes the microprocessor (IC16 pin 7) and the display module to be reset. The monitor then performs its powerup self-test routines, and if the “glitch” has been cleared, the monitor resumes normal operation. If the problem still exists, a self-test or other error should be displayed.

Serial I/O Controller

9.3.22

Digital data from the three Analog-To-Digital Convertors is read by the CPU through its clocked serial data input (RXS) at IC16 pin 52. The PEEL IC18 acts as the Clocked Serial I/O (CSI/O) Controller.

Except during powerup or Watchdog Timer reset, IC39 pin 3 provides an interrupt to the CSI/O controller in the form of a 5 millisecond period square-wave input to IC18 pin 7 (INT5MS).

On the rising edge of INT5MS, a CPU interrupt request is generated when IC18 pin 18 (CPUINIT*) goes low. The CPU responds by sending the clock input to CSI/O controller (CKS) at IC18 pin 6 low. (This CKS line is inactive high unless a serial receive operation is in progress.) The CPU also sets up the ADC decode lines AA1 and AA0 at IC18 pins 5 and 4, and as a result, one of the ADC chip select lines (CSADC1*, CSADC2*, CSADC3*) is brought low, and the CPUINIT* line is disabled.

On the rising CKS signal a CLKS output pulse at IC18 pin 14 is sent as a serial clock input to the ADC selected by the decode lines. Decode results are shown below.

AA1	AA0	Decode
0	0	Red LED 20-bit ADC
0	1	Infrared LED 20-bit ADC
1	1	Sensor Status 8-bit ADC
1	0	Internal CSI/O signal (TEND)

Table 2. CSI/O Decode Lines

Successive CKS/CLKS pulses cause the ADC data to be shifted out of the ADC (most significant bit first) along the serial data line (SDATA) to the CPU serial input (RXS) at IC16 pin 52.

After receiving the correct number of bits for the ADC being read, the CPU changes the AA1 and AA0 decode lines and exerts the Next line (NEXT*) at IC18 pin 9 low. This restarts the serial data shifting out of the newly selected ADC.

After all three ADCs have been read, the CPU sets the AA1 and AA0 decode lines to exert the internal TEND signal and set the 8-bit ADC to the next channel (so that it has time to settle before the next read of the ADC). This re-enables the CPUINIT line. At this point the CSI/O controller is reset awaiting an INT5MS pulse to begin the cycle again.

RS232 Serial Communication

9.3.23

The Model 520A supports serial (RS232) communication with external devices via the monitor's rear panel connector. Signals to and from the rear panel RS232 connector are electrically isolated from the rest of the Model 520A electronics by four opto-isolators (IC3-IC6). (See sheet 5 on 2471 schematic, sheet 4 on 2775 schematic.)

An isolated secondary coil from transformer T1 is rectified and filtered by D1 and C1 before being input to the +5 volt regulator IC2. The regulated output of IC2 is sent to pin 25 of the 25-pin D connector on the rear panel in order to power the optional Analog Module or Telemetry transmitter. The supply also powers IC1.

The Dual RS232 Transmitter/Receiver, IC1, uses a single +5 volt supply (pin 16). The two Receive (Rx) inputs can accept ± 30 volt levels, while the two Transmit (Tx) lines output ± 9 volt levels. The four level translators within the chip turn the RS232 level signals to 5 volt TTL/CMOS compatible levels.

Two signals lines TX0 and TX1 transmit data from the CPU across the opto-isolators ICs 3 and 4 to IC1. Here the signals are level shifted to the standard ± 9 volt levels and sent to the rear panel connector.

When the CPU brings the TX0 signal at R62 Low, the LED portion of IC3 becomes forward biased and illuminates and causes the transistor portion of the device to conduct. Therefore, when TX0 is Low, IC1 pin 11 is Low and when TX0 is High, IC1 pin 11 is also High. All four opto-isolators work in the same way. IC1 level shifts its input at pin 11 and outputs RS232 standard ± 9 volt levels at the TX signal at J101 pin 3. This transmit signal is used to output data to RS232 devices including the Hewlett-Packard ThinkJet Printer and the Novamatrix Model 315 Printer.

The CPU TX1 signal crosses the isolation barrier at IC4. The 0-5 volt opto-isolator output at IC4 pin 5 is used to bring data to the Optional 9622-01 Analog Module. The IC4 pin 5 output is also brought to IC1 pin 10. The level shifted output at IC1 pin 7 is unused (except for factory test use).

The Receive (RX) line at J101 pin 2 and the Clear To Send (CTS) line at J101 pin 6 are input signals to the Model 520A. They are level shifted by IC1 and sent across the isolation barrier by ICs5 and 6 respectively.

The transmit signal TX1, is dedicated to communication with the Optional Analog Module (Catalog Number 9622-01) which when connected to the rear panel connector, provides analog representations of the SpO₂ and Pulse Rate values, a plethysmogram signal, and a pass through port for the RS232 connector.

The transmit output TX0 from the CPU and the Receive (RX0) and Clear To Send (CTS*) inputs to the CPU are connected to the rear panel RS232 connector.

10

Maintenance

General

10.1

This section presents recommended maintenance schedules for the Model 520A and information on general maintenance, such as battery and fuse replacement, disassembly and assembly instructions, and system software updates.

Maintenance Schedules

10.2

The electronic circuits within the Novamatrix Model 520A Pulse Oximeter monitor do not require scheduled calibration or service. However, in order to maximize battery life, the monitor's internal battery should be exercised monthly. Novamatrix recommends the following maintenance schedules.¹

- **Cleaning and Sterilization:**
Perform as required. See *Cleaning and Sterilization* on page 36.
- **Battery Life and Maintenance:**
See *Battery Life and Maintenance* on page 37.
- **Functional Test:**
The test may also be used as a “spot check” to verify system operation if reports of malfunctions are received. The test verifies overall functional integrity of the monitor and sensors. Section 12, *Functional Test*, on page 51
- **Accuracy Test:**
This test, which requires the use of the Model TB500B Sensor Simulator, verifies the performance accuracy of the Model 520A. If the monitor does not pass the accuracy test, the Calibration and Adjustment test should be performed. See Section 13, *Accuracy Test*, on page 55.
- **Calibration Tests:**
These tests contain information on calibrating the electronic circuits within the Model 520A and should only be performed if the monitor fails to pass the Functional and/or Accuracy Tests. Only qualified service personnel should attempt to perform the Calibration and Adjustment Test. See Section 14, *Calibration Tests*, on page 59.

1. At the customer's request, Novamatrix will provide repair and calibration services under the terms of a Service Contract. Contact the Novamatrix Service Department for contract details.

Cleaning and Sterilization

10.3

Model 520A Monitor

10.3.1

- Turn the monitor off and unplug from the AC (Mains) line before cleaning.
- The monitor can be cleaned and disinfected with solutions such as a 70% isopropyl alcohol, 2% gluteraldehyde, or 10% bleach solution. Then wipe down with a water-dampened clean cloth to rinse. Dry before use.
- Do not immerse the monitor.
- Do not attempt to sterilize the monitor.

Finger Sensor

10.3.2

- The sensor can be cleaned and disinfected with solutions such as a 70% isopropyl alcohol, 2% gluteraldehyde, or 10% bleach solution. Then wipe down with a water-dampened clean cloth to rinse. Dry before use.
- Make certain that the finger sensor windows are clean and dry before reuse.
- Do not immerse the finger sensor.
- Do not attempt to sterilize the finger sensor.
- After cleaning the finger sensor, verify that the sensor is physically intact, with no broken or frayed wires or damaged parts. Make certain that the connectors are clean and dry, with no signs of contamination or corrosion. Do not use a broken or damaged sensor or one with wet, contaminated or corroded connectors.

Y-SENSOR™ and Y-STRIP™ Taping System

10.3.3

- Do not immerse connector on the Y-Sensor.
- The Y-Sensor may be immersed—up to, but not including, the connector, in a 2% gluteraldehyde solution, or 10% bleach solution. Refer to manufacturer's instructions and standard hospital protocols to determine recommended times for disinfection and sterilization.
- Rinse thoroughly with water and dry before use (do not rinse the connector).
- Do not attempt to sterilize Y-Sensor except as stated above.
- After cleaning or sterilizing the Y-Sensor, verify that the sensor is physically intact, with no broken or frayed wires or damaged parts. Make certain that the connectors are clean and dry, with no signs of contamination or corrosion. Do not use a broken or damaged sensor or one with wet, contaminated, or corroded connectors.
- Treat Y-Strip Tapes and foam wraps in accordance with hospital protocol for single-patient use items.

Battery Life and Maintenance

10.4

The monitor was shipped from the factory with a fully charged battery. Since the monitor draws some battery power even while turned off, it is recommended that the monitor be operated on line power for the first 24 hours to allow ample time for the battery to fully recharge.

The monitor may be operated on line power while the battery is recharging. Approximately 16 hours are required to fully recharge a discharged battery.

During prolonged periods of storage or shipment, the battery may discharge enough to prevent the unit from turning on while operated from battery power. If this occurs, plug in the line cord, set the switch on the rear panel to the “I” On position, ensure that the green \sim indicator on the front panel is illuminated, and allow the monitor to charge for 24 hours before switching it on.

Mains Voltage Configuration

10.5

The rear panel power entry module indicates the mains voltage setting for the monitor. Check that the voltage is correct before attaching the AC line cord and powering the monitor. The Model 2001 can be set to operate from 100-120 VAC 50/60Hz or 200-240VAC 50/60Hz.

Instructions for fuse replacement and changing the mains voltage setting follow.


CAUTION: Replace fuses with same type and rating. Verify proper fuse value for mains voltage setting (see table below).

Fuse Replacement


10.5.1

1. Check that the monitor is OFF.
2. Set the rear panel power entry module switch to OFF (“O”). Remove the AC line cord from the power entry module.
3. Using a flat blade screwdriver, pry the fuse access door open to expose the fuse housing. Note the orientation of the fuse housing (this determines the mains operating voltage).


Mains Voltage	Fuses (Slo Blo)
100-120 VAC	0.5 A 250V
200-240 VAC	250mA 250V



4. Pry the fuse housing out from the power entry module.



5. Replace the blown fuse(s) with the proper type and rating.
6. Reinstall the fuse housing. When positioning the housing into the power entry module make sure that it is oriented correctly. Press the fuse housing back into the power entry module.
7. Close the fuse access door and verify that the proper mains operating voltage is displayed.



Changing the Mains Voltage Setting

10.5.2


1. Check that the monitor is OFF.
2. Set the rear panel power entry module switch to OFF (“O”). Remove the line cord from the power entry module.

CAUTION: Replace fuses with same type and rating. Verify proper fuse value for mains voltage setting (see table below).


Mains Voltage	Fuses (Slo Blo)
100-120 VAC	0.5 A 250V

200-240 VAC	250mA 250V
-------------	------------


- Using a flat blade screwdriver, pry the fuse access door open to expose the fuse housing. Pry the fuse housing out from the power entry module.



- Install the proper type and rating fuse for the mains voltage setting required.
- Position the housing into the power entry module so that the desired voltage is furthest away from the switch (see below).



- Close the fuse access door and verify that the proper mains operating voltage is displayed.




AC Mains on older style assemblies

10.5.3


The Model 520A will operate on line voltages of 100, 120, 220 or 240 VAC \pm 10% at 50/60 Hz. At 100/120 VAC, two 0.5 Amp, 250 Volt, Slo-Blo (time delay) AC Mains fuses are required, while Two T 250 mA, 250 Volt fuses are required when operating at 220 or 240 VAC. For continued protection against fire hazard, replace only with fuse(s) of the same type and rating. The rear panel Power Entry Module (PEM) houses the AC Mains fuses and is where the input voltage is selected. It also houses the AC Mains Power switch: set to “I”, AC Mains power is presented to the internal circuitry; set to the “0” position, power does not pass through the PEM.

Replacing the AC Mains Fuse(s)**10.5.4**

1. Turn off the Model 520A. Set the rear panel AC Mains power switch to “0” and disconnect the line cord from the monitor.
2. Place a screwdriver into the Fuse Access Slot and pry open the Fuse Access Door.

**Figure 4.** Power Entry Fuse Access Door Opening

3. With the Fuse Access Door open pull the fuse(s) out by pulling on the point of the arrow indicator, the fuse holder assembly will slide out.

**Figure 5.** Fuse Removal

4. After replacing the fuse(s), snap the fuseholder assemblies into the PEM and shut the Fuse Access Door.

Changing the AC Mains Voltage**10.5.5**

1. Turn off the Model 520A. Set the rear panel AC Mains power switch to “0” and disconnect the line cord from the monitor.
2. Place a screwdriver into the Fuse Access Slot and pry open the Fuse Access Door.

- Using needle-nosed pliers, pull the Voltage Selector Drum from the PEM. Note the orientation of the drum; the proper voltage should face out.




Figure 6. Power Entry Module Voltage Selector Drum Removal

- Set the Voltage Selector Drum so that the printed side of the correct voltage faces you. The voltage selections are pictured below.




Figure 7. Power Entry Module Voltage Selection Adjustment

- Snap the Voltage Selector Drum back into the PEM. Be sure that the Voltage Setting Indicator shows the proper voltage. Close the Fuse Access Door.

Assembly Exchanges

10.6

Disassembly should be performed by qualified personnel. Follow proper grounding procedures to avoid damage to internal components from static discharge.

1. Turn the Model 520A Off. Disconnect the line cord and sensor. Remove the four cover screws from the bottom cover. Holding both case halves together, flip the monitor right-side up.




Figure 8. Removing the Top Cover

2. Lift the top cover from the monitor. Use a gentle rocking motion to lift first one side and then the other side a little at a time. Set the Red Alert Bar lens aside with the cover for safe keeping.




Figure 9. Monitor Assembly

3. Unplug P102 The Power Supply connector from the 2471 Main Board, disconnect P109 (Speaker), P108 (Sensor), P104 (Keypanel), P106 (Backlight), P105 (Alert Bar), and P103 (Display). Unscrew 2471 Main

Board from the rear panel. Remove the front Bezel assembly then slide the 2471 Main Board out, the 2472 Power Supply Board and Battery should be exposed.

4. The separate sections of the monitor can now be removed.




Figure 10. 520A Assembly with 2471 Main Board Removed


5. **IMPORTANT.** The battery is connected to the Power Supply Board through connector J302. Before attempting to apply power to the monitor ensure all connections are properly made, then connect the AC mains and turn the PEM switch to “I”, ensure that the ~ on the front panel illuminates before powering up.
6. Remove/replace the various assemblies as needed. **IMPORTANT:** The monitor *will be damaged* if power is applied to it while cables or assemblies are improperly connected.
7. After the monitor is fully assembled perform safety checks. Reconnect the line cord to the AC Mains and to the monitor. Set the AC Mains switch to the “I” (On) position. Verify the ~ indicator illuminates. Turn the monitor on.
8. Measure the AC leakage current from the monitor’s chassis to earth ground with the monitor grounded, ungrounded, and ungrounded reverse polarity. When operating from 100/120 VAC, no leakage current may be greater than 25 μ A. If operating at 220/240 VAC, no leakage current may be greater than 50 μ A.
9. With the monitor grounded, measure from the AC line to the SaO₂ Input connector pins. When operating from 100/120 VAC, no leakage current may be greater than 25 μ A. If operating at 220/240 VAC, no leakage current may be greater than 50 μ A.

Changing System Software

10.7

The system software is contained in EPROM IC17 on the 2775 Main Board. New software releases are made available from time-to-time. These new releases may add features or be maintenance upgrades. To install a new EPROM:

1. Follow the steps listed in *Assembly Exchanges* on page 41 to open the monitor. Be sure to disconnect J102 from the 2775 Main Board to ensure that there is no power on the board while changing the system software.
2. Use a PLCC extraction tool to carefully pry the EPROM IC17 from the socket.



3. Install the update EPROM into the socket. Align the EPROM so that IC17 pin-1 is inserted into pin 1 of the socket.

Changing System Software on 2471 Main Board

10.7.1

The system software is contained in EPROM IC17 on the 2471 Main Board. New software releases are made available from time-to-time. These new releases may add features or be maintenance upgrades. To install a new EPROM:

1. Follow the steps listed in *Assembly Exchanges* on page 41 to open the monitor. Be sure to disconnect J102 from the 2471 Main Board to ensure that there is no power on the board while changing the system software.

2. Use a small flat-blade screwdriver (or IC extraction tool) to pry the EPROM IC17 from the socket—be careful not to bend the pins.




Figure 11. Changing the System Software EPROM

3. Install the new EPROM into the socket. Align the EPROM so that IC17 pin-1 is inserted into pin 1 of the socket.
4. Reverse the above steps to reinstall the assembly into the monitor.




[This page intentionally blank]

11

Troubleshooting

Fault, alert and error messages that can appear on the Model 520A Message Center display are described below.

Paragraphs marked with a “” are intended for qualified service personnel only.

Message Center Display	Possible Explanation
AUDIO OFF DISABLED	Displayed if user tries to enable Audio Off mode (by pressing and holding the AUDIO key) while the “Allow Audio Off” portion of the Options Menu is set to “No”.
BAD SIGNAL TIMEOUT	Monitor not receiving valid signals from sensor. May be caused by excessive motion, cardiac arrhythmia or other situations leading to poor signal. Check patient status, reposition sensor.
BATTERY EXHAUSTED CONNECT LINE CORD	Monitor is running on battery power and the battery power has been depleted. Connect line cord to AC Mains power source and set the rear panel switch to “I”.  Monitor’s rear panel fuse has blown, monitor switched over to battery power and has depleted battery life.
CAN’T I.D. PROBE	Sensor is placed on too thick of a site. Reposition the sensor on a thinner (less opaque) section of tissue. A non-SuperBright™ sensor is connected, use only 87xx series sensors. Sensor is faulty. Remove sensor from use and contact qualified service personnel.
CONNECT SP02 PROBE	Sensor is disconnected from the monitor. Sensor is faulty. Remove sensor from use and contact qualified service personnel.  Both sensor LEDs have failed or the cable connections are open or shorted.
PROBE FAULTY	Sensor faulty. Remove sensor from use and contact qualified service personnel.  This error traps intermittent faults such as Probe Red LED or Infrared LED failures caused by frayed or broken wires in the sensor cable. Whereas PROBE RED or I.R LED FAILED is displayed for as long as the error lasts, Faulty Probe is displayed if those errors cease without the user replacing the faulty sensor.
EVENT MARKED	An event was successfully entered into trend memory.

Message Center Display	Possible Explanation
INCOMPATIBLE PROBE	<p>A non-SuperBright™ sensor is connected, use only 87xx series sensors.</p> <p>Sensor is faulty. Remove sensor from use and contact qualified service personnel.</p>
INSUF. LIGHT **	<p>Insufficient Light, where ** is time in seconds (after 99 seconds display shows "--"). Sensor placed on a site too thick (or opaque) for adequate light transmission. Reposition the sensor.</p> <p>A non-SuperBright™ sensor is connected, use only 87xx series sensors.</p>
LIGHT INTERF. **	<p>Light Interference, where ** is the time in seconds (after 99 seconds display shows "--") ambient light sources (sunlight, warming lights, etc.) are interfering with sensor light sources. Shield the sensor from ambient light sources.</p>
LOW SIGNAL **	<p>Low Signal Strength, where ** is the time in seconds (after 99 seconds display shows "--") that the pulse strength as detected by sensor is too weak for proper monitor operation. Reposition sensor.</p>
MONITOR ERROR	<p>Monitor faulty. Record error message (appearing on bottom line of display) and contact qualified service personnel.</p> <ul style="list-style-type: none"> ☞ "RAM SELF TEST FAILED" - Ram failed the power up self-test. ☞ "ROM SELF TEST FAILED" - Calculated EPROM checksum not does not equal stored value. Bad EPROM. ☞ "BAD STACK POINTER" - Stack pointer base not at top of stack. ☞ "CORRUPT SYSTEM MMU" - The CPU MMU pointing to an illegal address map. Digital Bd fault. ☞ "STACK OVERFLOW" - Stack pointer exceeded allocated stack size. ☞ "DISP BUFFER OVERFLOW" - Display buffer queue exceeded its allocated size. ☞ "HARDWARE ERROR 1" - (1)Main Bd PEEL is defective, (2) Main Bd 20-bit ADCs exceeded acceptable calibration thresholds. ☞ "XXms INT.TIMEOUT" where XX is 10 or 25. The interrupt has interrupted itself. ☞ "ERROR UNDEFINED" - Failed in an unknown state.
MONITOR PERFORMING SELF TEST.	<p>Monitor is performing its power up system diagnostic tests.</p>
Parameters Reset To Factory Default	<p>Displayed when monitor is turned on while pressing the ALERT RESET key. Monitor now using factory default settings.</p>

Message Center Display	Possible Explanation
PROBE FAULTY RD	Sensor faulty. Remove sensor from use and contact qualified service personnel. ☞ Sensor Red LED has failed or the cable connections are open or shorted.
PROBE FAULTY IR	Sensor faulty. Remove sensor from use and contact qualified service personnel. ☞ Sensor Red LED has failed or the cable connections are open or shorted.
PROBE OFF PATIENT	Sensor disconnected from patient, improperly applied, or placed on an area too translucent for proper sensor operation. Reposition sensor.
PULSE-HIGH	Selected pulse rate high alert limit has been violated.
PULSE-LOW	Selected pulse rate low alert limit has been violated.
PULSE OUT OF RANGE	Pulse rate is less than 30 bpm or is greater than 250 bpm.
Revision 2.0 APR/28/92	Monitor software revision level. Displayed when monitor is turned on while pressing AUDIO key.
SpO2-HIGH	Selected saturation high alert limit has been violated.
SpO2-LOW	Selected saturation low alert limit has been violated.

[This page intentionally blank]

12

Functional Test

Introduction

12.1

The test verifies overall functional integrity¹ of the monitor and sensors. If the monitor or sensors do not pass these tests, remove from use and contact the Novamatrix Service Department for repair/replacement assistance.

This procedure assumes the technician performs each step as indicated - leaving the monitor in a known state prior to performing the next step. If steps are omitted or performed out of order, be sure that the monitor is set to the correct state before continuing.

Monitor Functional Test

12.2

- 12.2.1 Equipment Required:
Model 520A to be tested
AC Mains Line Cord (supplied with monitor)
Finger Sensor (Cat. No. 8776) or Y-Sensor™ (Cat. No. 8791)
- 12.2.2 Visually inspect the monitor and line cord. Verify there is no external damage. Shake the unit to check that there is no loose hardware inside.
- 12.2.3 Inspect the sensor(s). Verify mechanical integrity.
- 12.2.4 Check the rear panel power entry module. Verify that it is set to the proper voltage.
- 12.2.5 Set the AC Mains switch to the “0” (Off) position. Connect the line cord to the monitor and to the AC Mains.
- 12.2.6 Verify the ~ (AC Mains) indicator is not illuminated.
- 12.2.7 Set the AC Mains switch to the “I” (On) position. Verify the ~ indicator illuminates.
- 12.2.8 Depress and hold the **ALERT RESET** key. Press the front panel **POWER** key then release both keys. The monitor will;
- Briefly illuminate all displays and indicators²
 - Produce an audible “beep” indicating the audio is operational




1. See Section 13, *Accuracy Test*, on page 55 for tests that check the accuracy of the monitor.

2. The AC indicator will illuminate only when the monitor is connected to line power and the rear panel Power switch is On (I).

- Briefly display Parameters Reset To Factory Default.³
 - Briefly display MONITOR PERFORMING SELF TEST.⁴
 - The message center will display ERASE STORED TRENDS ? for approximately 10 seconds.
- 12.2.9 Verify the Message Center flashes CONNECT SPO2 PROBE while no sensor is connected to the monitor.
- 12.2.10 Press the **ALERT** key to display the SET ALERT LIMITS menu. Press the ↓ (down arrow) key several times and verify the Saturation display upper alert limit decreases each time the key is pressed.
- 12.2.11 Press the **RUN** key. Note the displayed Saturation alert limit settings. Turn the monitor off and back on again. After the self-test is complete, verify the monitor “remembered” the alert limit values.
- 12.2.12 Press the **Menu** key and then the **AUDIO** key. Press the **ALERT** key to display the SET ALERT VOLUME menu. Press and hold the ↓ key and verify the volume of the tone decreases as the volume display decreases to 01. Press the ↑ key to return to maximum (07) volume. Press **RUN**.
- 12.2.13 Press **MENU** then press **LITE**, verify the backlight toggles between dim and bright with each depression. Press **RUN**.
- 12.2.14 Press the **Menu** key and then the **AUDIO** key. Press the **PULSE** key to display the SET AUDIO FEATURES menu. Press the ↑ key to set the value to 01 then press **RUN**.
- 12.2.15 Do not apply the sensor to yourself. If using a Y-Sensor, fold the sensor heads so that the windows face (or even touch) each other. Connect the sensor to the monitor. Verify the Message Center flashes PROBE OFF PATIENT and no alarm sounds.
- 12.2.16 If using a Finger sensor, apply it to your index finger, and if using a Y-Sensor, tape the sensor to your index finger. Verify that after several seconds, reasonable Pulse Rate and Saturation (>95%) values are displayed.
- 12.2.17 Check that the Message Center displays a Pulse Activity Bar and that a tone sounds with each pulse beat. (Note that the pitch of the tone will vary with the Saturation display value.)
- 12.2.18 Remove the sensor from your finger. Verify the Saturation and Pulse Rate displays blank out, an alarm sounds and PROBE OFF PATIENT is displayed. (If using a Y-Sensor, you may have to fold the heads so that the windows face each other in order to make the alarm and message activate.)
-

3. This message is only displayed if the ALERT RESET key is depressed at power up.

4. If “Monitor Error” appears, turn the monitor off and back on. If the message reappears, contact qualified service personnel. See Section 11, *Troubleshooting*, on page 47.

- 12.2.19 Press and release the **AUDIO** key. The  (Two Minute Silence) indicator should illuminate and the alarm should be silenced. After two minutes the indicator turns off and the alarm sounds again (providing the alarm condition still exists).
- 12.2.20 With the alarm sounding, press and hold the **AUDIO** key until the  (Audio Off) indicator illuminates. Verify that the alarm is silenced.
- 12.2.21 Press and release the **AUDIO** key. Verify the  indicator shuts off and the alarm sounds.
- 12.2.22 With the PROBE OFF PATIENT alarm still sounding, press the **ALERT RESET** key. Verify the message and alert indicators continue to flash but the alarm is silenced.
- 12.2.23 Set the AC Mains switch to the “0” (Off) position. Verify the unit operates from battery power (provided the battery is substantially charged). Place the AC Mains switch back to the “I” position.
- 12.2.24 Press the **POWER** key to turn the monitor off and disconnect the sensor.
- 12.2.25 Depress and hold the **ALERT RESET** key. Press the front panel **POWER** key then release both keys. The monitor will power up with the factory default settings in place. After the self-test is complete, turn the monitor off.
- 12.2.26 This completes the Functional Test. If the monitor and sensor performed as described above, they are functionally operational. If the monitor is to be returned to clinical use, be sure to let the user know that the monitor is now using its default settings—as these may differ from the user’s “normal” settings.

Special Power Up Functions

12.3

Each of the Model 520A front panel pushbutton keys is linked to a special “hidden” power up function. To initiate the special function, start with the monitor off and press a front panel key. While still holding that key, press the **POWER** key to turn the monitor on.

AUDIO = Software (EPROM) Revision Level

Turn the monitor on while holding the **AUDIO** key and the monitor’s software revision level is displayed in the message center for as long as the **AUDIO** is pressed. Release the key and the monitor continues with its normal power up sequence.

ALERT RESET = Return to factory default settings

Turn the monitor on while holding the **ALERT RESET** key and the monitor resets all its control settings (including alert limits, averaging times, etc.) to the factory default values. A Parameters Reset To Factory Default message is briefly displayed and the monitor continues with its

normal power up sequence. Note that the system date and time are not affected by this action.




Figure 12. Softkey Identification

Softkey #1 = EPROM Checksum Test, Production Test Mode

Turn the monitor on while holding **Softkey #1** (the leftmost softkey) and the monitor calculates and displays the checksum value of the EPROM containing the system software. Press and hold **Softkey #3** while the checksum is displayed and the monitor enters Production Test Mode.

Softkey #2 =Extended Display and Watchdog Test

Turn the monitor on while holding **Softkey #2** (second from left) and all front panel indicators and display segments illuminate for several seconds. The watchdog now resets the monitor which should then execute its normal power up sequence. Note that the ~ indicator will illuminate only if the monitor is being powered from the AC Mains.

13

Accuracy Test

Introduction

13.1

The Accuracy Test verifies the performance accuracy of the Model 520A Pulse Oximeter monitor. If the monitor does not pass the accuracy test, the monitor should be recalibrated. Refer to Section 14, *Calibration Tests*, on page 59, or contact the Novametrix Service Department for repair/replacement assistance.

The test requires the use of the Model TB500B Sensor Simulator.¹ This is the same device used by the factory technicians to calibrate the monitor prior to shipping. The TB500B is an updated version of the TB500A Test Box. Owners of TB500A should contact the Novametrix Service Department for details on upgrading to the TB500B. Note that the TB500A, used in conjunction with the Cat. No. 5453-00 Adapter Cable, may be substituted for the TB500B in most parts of this test.

Monitor Accuracy Test


13.2

- 13.2.1 Equipment Required:
Model 520A to be tested
AC Mains Line Cord (supplied with monitor)
TB500B Sensor Simulator (Cat. No. 5530-00)
- 13.2.2 Press and release the front panel **POWER** key. The monitor will;
- Briefly illuminate all displays and indicators²
 - Produce an audible “beep” indicating the audio is operational
 - Briefly display MONITOR PERFORMING SELF TEST.³
 - Display the Message Center flashes ERASE STORED TRENDS? for ten seconds then displays CONNECT SPO2 PROBE while no sensor is connected to the monitor.
- 13.2.3 Press the **Menu** key and then press **AVG**.
The SELECT SPO2 AVERAGING menu should be displayed. Press **2s** and then press **RUN** to return to the Main Menu. This sets the SpO₂ averaging to 2-seconds for the quickest response to changing TB500B settings.

1. Available through the Novametrix Service Department.

2. The AC indicator will illuminate only when the monitor is connected to line power and the rear panel Power switch is On (I).

3. If “Monitor Error” appears, turn the monitor off and back on. If the message reappears, contact qualified service personnel. Refer to Section 11, *Troubleshooting*, on page 47.

- 13.2.4 Set the TB500B to these settings;
Power - On, Sensor Type⁴ - 87xx, Signal Attenuation - 3, Saturation - 92.
- 13.2.5 Connect the TB500B to the Model 520A SpO₂ Input Connector.
Verify a Pulse Activity Bar appears in the Message Center.
Verify SpO₂ and Pulse Rate values appear after several “pulses”.
- 13.2.6 Set the Signal Attenuation to 1. Verify the Saturation and Pulse Rate displays blank out, an alarm sounds and PROBE OFF PATIENT is displayed. Set the Signal Attenuation to 3 and verify the displays and Pulse Activity Bar return.
- 13.2.7 Set the Saturation to 0. Verify that after a short delay, LOW SIGNAL appears in the Message Center. Set the Saturation to 100 and verify the displays and Pulse Activity Bar return.
- 13.2.8 Press and hold the TB500B RED push-button. Verify the monitor alarms and displays PROBE FAULTY RD. Release the button.
- 13.2.9 Press and hold the TB500B INFRARED push-button. Verify the monitor alarms and displays PROBE FAULTY IR. Release the button. Verify the displays and Pulse Activity Bar return.
- 13.2.10 Press and hold both the TB500B RED and INFRARED push-buttons. Verify the displays blank out and CONNECT SPO2 PROBE is displayed. Release the buttons. Verify the displays and Pulse Activity Bar return.
- 13.2.11 Disconnect the TB500B from the monitor. Set the Sensor Type switch to 86xx.⁵ Reconnect the simulator to the monitor. Verify the message INSUFFICIENT LIGHT is replaced in approximately 10 seconds by CAN'T I.D. PROBE, which is itself replaced later by PROBE FAULTY. Disconnect the TB500B from the monitor, return the Sensor Type switch to 87xx and then plug it back in. Verify the displays and Pulse Activity Bar return.
- 13.2.12 Press and hold the **AUDIO** key until the  (Audio Off) indicator illuminates. This will keep the monitor silent for the remainder of the test.
- 13.2.13 Verify the displayed SpO₂ value is within the tolerances stated below for each setting of the Signal Attenuation and Saturation switches. Verify a Pulse Rate of 60 bpm ± 1bpm for all switch settings. (Note that alert

4. The 5453-00 Adapter Cable is used instead of the Sensor Type switch for the TB500A.

5. If using a TB500A, disconnect the Adapter Cable and plug the TB500A directly into the Model 520A. The Message Center should display “INCOMPATIBLE PROBE”. Install the Adapter Cable between the TB500A and the Model 520A and move on to the next step.

messages will be generated and displayed as the saturation value violates the alert limit settings.)

TB500B		520A	520A SpO ₂
SAT	ATTEN	SpO ₂	(w/ TB500A)
100	3	99 ± 1	99 ± 1
92	3	92 ± 2	92 ± 2
82	3	82 ± 2	84 ± 2
72	3	72 ± 2	77 ± 2
62	3	62 ± 2	69 ± 2
100	7	99 ± 1	99 ± 1
92	7	92 ± 4	92 ± 4
82	7	82 ± 4	84 ± 4
72	7	72 ± 4	77 ± 4
62	7	N/A	N/A

Table 3. SpO₂ Display tolerances for TB500B settings

- 13.2.14 Press the **Menu** key and then press **AVG**. Press the appropriate key to return the averaging to its pretest value.
Press **RUN** to return to the Main Menu.
- 13.2.15 This completes the Monitor Accuracy Test. If the monitor does not meet the above listed specifications, refer to Section 14, *Calibration Tests*, on page 59, or contact the Novamatrix Service Department for recalibration, repair, or replacement information.

[This page intentionally blank]

14

Calibration Tests

Introduction

14.1

The Calibration Tests verify the calibration and operation of the electronic circuits within the Model 520A Pulse Oximeter. These tests do NOT need to be performed on a regular (preventative maintenance) basis. Perform these tests only if the monitor fails to operate as expected and/or fails the *Functional Test* on page 51 and/or the *Accuracy Test* on page 55. Calibration Tests should be performed only by qualified service personnel. The Model 520A contains static sensitive devices. Follow proper grounding procedures when handling the internal components to avoid damage from static discharge.

The test requires the use of the Model TB500B¹ Sensor Simulator Test Box. This is the same device used by the factory technicians to calibrate the monitor prior to shipping. The TB500B is an updated version of the TB500A² Test Box.

If the monitor does not pass the Calibration Tests, remove it from use and contact the Novamatrix Service Department for repair/replacement assistance.

This procedure assumes the technician performs each step as indicated—leaving the monitor in a known state prior to performing the next step. If steps are omitted or performed out of order, be sure that the monitor is set to the correct state before continuing.

Equipment Required and Test Setup

14.2

- 14.2.1 Model 520A to be tested
AC Mains Line Cord (supplied with monitor)
TB500B Sensor Simulator (Cat. No. 5530-00)
or TB500A Sensor Simulator and Adapter Cable 5453-00.
RS232 test fixture, PN: 5479-01
Saturation test jack, PN: 5942-00
Digital Voltmeter
Small Flat-blade screwdriver
- 14.2.2 With the Model 520A Off. Disconnect any line cord and sensor from the monitor.

1. Available through the Novamatrix Service Department.

2. The TB500A, used in conjunction with the Cat. No. 5453-00 Adapter Cable, may be substituted for the TB500B in most parts of this test.

- 14.2.3 Flip the monitor over to expose the bottom cover and remove the four cover screws. Holding both the top and bottom halves together, flip the monitor over again and set it on its bottom cover.
- 14.2.4 Lift the top cover from the monitor and set it aside. Set the Red Alert Bar lens aside with the cover for safe keeping.

Test Procedure


14.3

Unless otherwise specified all measurements are with respect to isolated ground. Use the top of R4 (end closest to IC2) as ground reference.

- 14.3.1 Attach a line cord and place the mains switch to the “I” position. Verify that the green \sim LED on the front panel illuminates. The monitor may power up when the Mains switch is set to “I”, if this occurs turn the monitor off with the front panel button.
- 14.3.2 Measure the voltage at J102 pin 2 and verify $13.2 \text{ VDC} \pm 0.4 \text{ volts}$. This voltage is dependent on a fully charged battery installed. If an installed battery is less than fully charged, the measured voltage will not match the expected voltage.
- 14.3.3 Simultaneously press the leftmost softkey and the power button. While Calculating Checksum is displayed, press and hold the third softkey from the right until Novamatrix Inc Production Test appears.
- 14.3.4 Measure the voltage at TP1. Verify $5.000\text{V} \pm 50\text{mV}$. Adjust VR1 if necessary.
- 14.3.5 Verify the following voltages:

<u>Location</u>	<u>Signal Name</u>	<u>Voltage</u>
IC40 pin 4	-VA	$-12.00\text{V} \pm 500\text{mV}$
IC40 pin 8	+VA	$12.00\text{V} \pm 500\text{mV}$
IC20 pin 28	VBACK	$4.75\text{V} \pm 250\text{mV}$
IC37 pin 14	+V5	$5.00\text{V} \pm 250\text{mV}$
IC37 pin 7	-V5	$-5.00\text{V} \pm 250\text{mV}$
IC37 pin 10	VREF	$2.500\text{V} \pm 5\text{mV}$
IC36 pin 1	4VREF	$4.045\text{V} \pm 50\text{mV}$
IC30 pin 14	-4VREF	$-4.045\text{V} \pm 100\text{mV}$

- 14.3.6 Short out capacitor C90 by installing a jumper wire from R59 to R56.



- 14.3.7 Connect the RS232 test fixture to J101 and the saturation test jack to the saturation input connector.
- 14.3.8 Run monitor test 1. Press the Test softkey, Automatic Test will appear. Press the Test softkey again to initiate the test, follow the prompts in each test and verify each test passes.

*Note: A series of auto tests will be performed, if any failures are encountered the monitor will default to MONITOR TEST 2. Of the list below only the tests with * next to them are performed in the auto test. Upon completion of the auto test the monitor will default to MONITOR TEST 2. The following is a list of the monitor tests available;*

Test#	Description
1	Auto Test
2*	RAM Test
3*	ROM Test
4*	Real Time Clock Test
5*	Audio Volume Test
6*	Audio Frequency Test
7*	RS232 Loopback Test
8	Keyboard Test
9	Mains Test
10*	20 Bit ADC Test
11*	8 Bit ADC Test
12	Display ADC Channels Test
13	515A-> PC Interface Test

- 14.3.9 Disconnect the jumper wire, the RS232 test fixture, and the saturation test jack.
- 14.3.10 Connect the TB500B to the saturation input connector. Set the saturation to 0 and the attenuation to 5.

14.3.11 Using the Next softkey increment to monitor test 12 (Display ADC channels). Press the Test softkey, Display ADC channels will appear. Press the Test softkey to initiate the test. Press the Prev softkey. Verify the LED pwr=20. Press the Prev softkey until IR Adc=xxxxxx appears, verify the value displayed is within the limits specified below. Press the Prev softkey again, RED Adc=xxxxxx will appear. Verify the value displayed is within the limits specified. Press the exit softkey and set the TB500B saturation to 100 and attenuation to 3.

RED ADC	5B32-88CB hex
IR ADC	3324-4CB7 hex

14.3.12 Power the monitor down.

14.3.13 Power up the monitor by pressing the power button. Verify a Monitor performing self test message is displayed then is immediately followed by the main monitoring screen.

14.3.14 Press the Menu softkey followed by the Avg and 2s softkey then press Run.

14.3.15 Set the TB500B to the settings in the chart. Verify the saturation values are within the specified tolerance:

TB500B Settings		Board Parameters	
<u>Saturation</u>	<u>Attenuation</u>	<u>Saturation</u>	<u>Pulse Rate</u>
100	3	98 - 100	60 ± 1
82	3	80 - 84	60 ± 1
62	3	60 - 64	60 ± 1
72	7	68 - 76	60 ± 1
92	7	88 - 96	60 ± 1

Note: Verify a visual (alert bar flashing) and audible alert condition occurs when the saturation is set to 82. Press and hold the audio alert key. Verify the audio off LED turns on and the audible alert is silenced.

14.3.16 Change the TB500B power switch to OFF. Verify an Insufficient Light error message is displayed.

14.3.17 Change the TB500B power switch to ON. Verify the error message clears and a saturation and pulse value are displayed.

14.3.18 Set the TB500B Signal Attenuation to “1”. Verify a Probe Off Patient message is displayed.

14.3.19 Set the TB500B Signal Attenuation to “3”. Verify the error message clears and a saturation and pulse value are displayed.

- 14.3.20 Set the TB500B Saturation to “0”. Verify a Low Signal Strength message is displayed.
- 14.3.21 Set the TB500B Saturation to “100”. Verify the error message clears and a saturation and pulse value are displayed.
- 14.3.22 Press the TB500B Red Open button. Verify a Probe Faulty message is displayed.
- 14.3.23 Release the Red Open button. Verify the error message clears and a saturation and pulse value are displayed.
- 14.3.24 Press the TB500B Infrared Open button. Verify a Probe Faulty message is displayed.
- 14.3.25 Release the Infrared Open button. Verify the error message clears and a saturation and pulse value are displayed.
- 14.3.26 Disconnect the TB500B from the saturation input connector. Verify a Connect SPO2 Probe message is displayed.
- 14.3.27 Set the current time / date. Press and hold the Menu softkey. Press the Next softkey until Monitor Options 2 appears. Press the Clock softkey. Use the appropriate softkeys to set the correct Time, Mins, Day, Month and Year.
- 14.3.28 Set the mains switch from “I” to “O”, verify the green ~ LED on the front panel goes out. Verify that the monitor continues to function on battery power without interruption. Power down the monitor by pressing the power button.

Safety Testing

14.4

- 14.4.1 Apply 1.5 KV for 60 seconds between the shorted hot and neutral leads of the power cord and chassis ground.
- 14.4.2 Apply 4 KV for 60 seconds between the shorted hot and neutral leads of the power cord and shorted saturation test jack.
- 14.4.3 Apply 1.5 KV for 60 seconds between the shorted saturation test jack and the chassis.
- 14.4.4 Measure the leakage from chassis gnd. to ref. gnd. with the monitor;
 -grounded
 -ungrounded
 -polarities reversed ungrounded
 With the monitor grounded measure the leakage from shorted probe to 120VAC 60Hz

[This page intentionally blank]

15

Connecting to other Equipment

The Model 520A Pulse Oximeter communicates with external equipment via its rear panel (25 pin D style female) connector. The monitor can be connected to RS232 compatible (digital) equipment. If the optional Analog Output Module is connected, the monitor can output to analog devices such as strip chart recorders. Refer to the Model 520A Pulse Oximeter User's Manual (Catalog No. 5693-23) for further information on the RS232 and analog outputs listed below.

Connecting the Seiko DPU-414 Thermal Printer 15.1

To connect a Seiko DPU-414 Thermal Printer (PN: 9140-00) to the *OXYPLETH*:

1. Connect the 9085-00 interface cable to the monitor's RS232 connector and to the printer's serial input connector—the 9 pin D connector.
2. Connect the AC adapter and turn the printer on.

Configuring the Seiko DPU-414 Printer 15.1.1

The Seiko DPU-414 Thermal Printer (Cat. No. 9140-00) must be configured to communicate with the *OXYPLETH*. When properly configured, the Seiko printer will retain the settings, even when turned off.

1. Slide the printer's power switch to OFF "O".
2. Press and hold the **ON LINE** button, then slide the power switch ON "I". Release the **ON LINE** button after the list of current settings starts printing out.

Setting the DIP switches: 15.1.2

1. The printout of the current settings is followed by the prompt:

"Continue? : Push On-line SW"

"Write? : Push Paper feed SW"

To change the DIP switch settings, push the **ON LINE** button (to leave the DIP switch settings unchanged, push the **FEED** button).

2. “DIP SW-1” will print. Enter the new settings for switches 1-8.

“ON” is set by pushing the **ON LINE** button once

“OFF” is set by pushing the **FEED** button once

The printer will confirm each selection. Repeat for DIP SW 2 and 3.

The DIP switch settings for the Model 520A are as follows:

```
[ DIP SW settings mode ]

Dip SW-1
  1 (OFF) : Input = Serial
  2 (ON ) : Printing Speed = High
  3 (ON ) : Auto Loading = ON
  4 (OFF) : Auto LF = OFF
  5 (ON ) : Setting Command = Enable
  6 (OFF) : Printing
  7 (ON ) : Density
  8 (ON ) : = 100 %

Dip SW-2
  1 (ON ) : Printing Columns = 40
  2 (ON ) : User Font Back-up = ON
  3 (ON ) : Character Select = Normal
  4 (ON ) : Zero = Normal
  5 (ON ) : International
  6 (ON ) : Character
  7 (ON ) : Set
  8 (OFF) : = U.S.A.

Dip SW-3
  1 (ON ) : Data Length = 8 bits
  2 (ON ) : Data Parity = No
  3 (ON ) : Parity Condition = Odd
  4 (ON ) : Busy Control = H/W Busy
  5 (OFF) : Baud
  6 (ON ) : Rate
  7 (ON ) : Select
  8 (ON ) : = 9600 bps

Continue ? : Push 'On-line SW'
Write ?   : Push 'Paper feed SW'
```

CAUTION: DIP SW Set Mode cannot be cancelled once it is initiated. Answer “ON” or “OFF” for every setting.

Note: More information about DIP switch settings can be found in the Seiko “DPU-414 Thermal Printer Operation Manual.”

3. When the printer finishes writing the new settings to memory, “DIP SW setting complete!!” is printed out and the printer returns to ON LINE mode.

CAUTION: Never turn the printer off while it is writing the new settings to memory. Wait until “DIP SW setting complete!!” is printed, then the printer power may be turned off.

Connecting the ThinkJet Printer

15.2

To connect the Hewlett-Packard ThinkJet Printer to the Model 520A:

1. The Hewlett-Packard ThinkJet must be a Model 2225D (RS-232C interface).
2. Set the dip switches on the rear panel of the ThinkJet.
Mode switches = 1, 2, 5, 6 up (on) and 3, 4, 7, 8 down (off).
RS-232C switches = 1 up (on) and 2, 3, 4, 5 down (off).




Figure 13. Connecting the ThinkJet Printer

3. Connect the interface cable, PN:5331-00, to the Model 520A RS232C connector and to the ThinkJet. (If Model 520A has an Analog Module attached, connect to the module's RS232C connector.)
4. Connect the printer's AC input and turn the printer on.
5. Select **T-Jet** (ThinkJet) from the **Ser.** (Serial) portion of the Model 520A Options Menu.
6. Refer to the ThinkJet Owner's Manual for further printer details.

Connecting the Model 315 Printer

15.3

To connect the Model 315 Printer to the Model 520A:

1. Set the dip switches on the Model 315 rear panel.
The settings are; 1-5 up (off), 6-7 down (on), and 8 up (off).

2. Connect the printer interface cable, PN:4913-00, to the RS232C connector¹ on the Model 520A. Connect the other end of the cable to the Model 315.




Figure 14. Connecting the Model 315 Printer

3. Connect the printer's AC adapter (if desired).
4. *Turn the printer on first*, then turn the Model 520A on.
5. Select **315** from the **Ser.** (Serial) portion of the Model 520A Options Menu.
6. Refer to the Model 315 documentation for further printer details.

Connecting Seiko DPU-411 Thermal Printer

15.4

To connect the SIEKO Model DPU-411 Thermal Printer to the Model 520A:

1. Set the DIP switches located on the bottom of the DPU-411.
 First bank of eight switches; 1,2 down (off), 3,4 up (on), 5,6 down (off), 7,8 up (on).
 Second bank of six switches; 1-3 up (on), 4-6 down (off).




Figure 15. Configuration of Seiko Thermal Printer.

2. Connect the printer interface cable, PN:5861-00, to the RS232C connector on the Model 520A. Connect the 25 pin connector at the other end to the 25 pin D connector on the rear panel of the DPU-411 (upper connector). Plug the cable DC power plug into the rear jack of the Model DPU-411.




Figure 16. Connecting the Seiko Model DPU-411 Printer.

3. Turn the printer on first, then turn the Model 520A on.
4. Select **Seiko** from the **Ser.** (Serial) portion of the Model 520A Monitor Options 2 Menu.
5. Refer to the Seiko Model DPU-411 documentation for further printer details.

Analog Output Module

15.5

The *optional* Analog Output Module (Catalog Number 9622-01), provides the necessary analog output voltages for use with analog instruments such as strip chart recorders. (The Model 520A does not directly support analog devices.)

The Analog Output Module attaches to the connector on the rear panel of the Model 520A and provides analog output voltages, via a 9-pin connector, and RS232C pass through, via a 25-pin connector, so that both analog and serial devices can be used simultaneously.

Two screws, supplied with the module, are passed through the module and screw into the pre-tapped holes in the rear panel of the Model 520A to secure the Analog




Figure 17. Analog Output Module

Output Module in place. The pinouts of the 9-pin analog and 25-pin RS232C connector are shown below.




Figure 18. Analog Output Module Pinouts


Analog Output Setup

15.6

A CALIBRATE RECORDER menu within the Model 520A software allows the user to easily calibrate analog recorders to the voltage levels produced by the Analog Output Module. To access this feature:

1. With the monitor in normal running mode press and hold the **MENU** key until SPO2 SETUP OPTIONS appears on the display. Press the **NEXT** key until


MONITOR OPTIONS 2 appears on the display. Select **REC** for CALIBRATE RECORDER menu.



Press and hold for 5 seconds while in normal running mode

Figure 19. Analog Output Setup

2. The Recorder Outputs menu will be displayed.



- **RUN** - return to Monitoring Mode
- **ZERO** - set analog outputs to 0 volts
(0% = 0 volts, 0 bpm = 0 volts, pleth = 0 volts)
- **HALF** - set analog output voltages to half scale
(50% = 0.50 volts, 125 bpm = 0.50 volts, pleth = 0.50 volts)
- **FULL** - set analog output voltages to full scale
(100% = 1.00 volt, 250 bpm = 1.00 volts, pleth = 1.00 volts)
- **PREV** - returns to MONITOR OPTIONS 2 menu.

3. Press the **ZERO** key and adjust the analog baselines on the recorder. Press the **FULL** key and adjust the recorder full scale deflection. Press the **HALF** key and check that the recorder registers at half scale.

Signal	Analog Connector Pin #	Min-Max Voltages (Referenced to Pins 2 + 3)		
		ZERO	HALF	FULL
Pulse Rate	1	0.000-0.012	0.477-0.520	0.955-1.030
Saturation	4 + 5	0.000-0.012	0.477-0.520	0.955-1.030
Plethysmogram	8 + 9	0.000-0.012	0.477-0.520	0.955-1.030

Table 4. Analog Output Module Voltage Tolerances

Note: Earlier versions of the Analog Output Module have a 2.5 volt maximum scale for the Plethysmogram and Pulse Rate (1.25 volt half scale).

4. When the recorder is properly adjusted, press the **RUN** key. The monitor will return to Monitoring Mode and automatically begin outputting analog data to the Analog Output Module.

Note: This setup procedure does not have to be performed each time the monitor is turned on—the monitor will always automatically output analog information via the Analog Output Module, this procedure simply helps to calibrate the recorder to the Model 520A signals.

16

Specifications

General

16.1

Specifications for the Novamatrix *OXYPLETH* Pulse Oximeter, Model 520A, are listed for informational purposes only, and are subject to change without notice.

Oxygen Saturation (SpO₂) Section

16.2

- Range, 0-100%
- Accuracy, 80-100% \pm 2%, 0-79% unspecified
(Approximately 68% of the observations are within the accuracy claim.)
- Display Resolution, 1%
- Averaging Time, menu selectable times of 2 and 8 seconds
- Audible SpO₂ Trend Feature
Pitch of (user selectable) Pulse Rate “beep” tracks the SpO₂ value (i.e., decreasing SpO₂ values are signalled by lower pitched “beeps”).
- Settling Time
Display settles to within 1% of the final reading less than 15 seconds after the sensor is properly applied.
- Alerts
Continuously displayed. Menu selectable high and low limits (100-50). Visible alarm is immediate. Audible alarm occurs after 10 seconds of continuous violation of the set limit, or immediately. Limit values are retained in memory when monitor is turned off, or the monitor can be set to use its default settings each time it is turned on.

Pulse Rate Section


16.3

- Range, 30-250 beats per minute (bpm)
- Accuracy, \pm 1% of full scale
(Approximately 68% of the observations are within the accuracy claim.)
- Display Resolution, 1 bpm
- Averaging Time, fixed at 8 seconds
- Settling Time
Display settles to within 1% of the final reading less than 15 seconds after the sensor is properly applied.

- Alerts
Continuously displayed. Menu selectable high and low limits (249-30 or Off). Visible alarm is immediate. Audible alarm occurs after 10 seconds of continuous violation of the set limit or immediately. Limit values are retained in memory when monitor is turned off, or the monitor can be set to use its default settings each time it is turned on.



General Specifications

16.4

- Operating Environment
50-104° F (10-40 °C), 0-90% relative humidity (non-condensing)
- Weight, 7 lbs 5 oz. (3.32 kg)
- Dimensions
Height, 3.3 inches (8.38 cm) Width, 9 inches (22.86 cm) Depth, 8 inches (20.32 cm)
- Power, 100-120/200-240 VAC, 50/60 Hz
- Fuse Rating
U.S.A.: 0.5 A, 250 V, Slo-Blo (x2) European: T 250 mA/250 V (x2)
- Battery
Type, lead-acid gel-cell Battery Life, 3 hours
(Note: Excessive alerting reduces battery life.) When 15 minutes of battery life remain, the  (low battery) indicator illuminates. When the battery becomes exhausted, the monitor display shuts down. Connect to AC power to recharge battery. Recharge Time, battery fully recharged in 12-15 hours max.

Additional Features

16.5

- 2 Minute Silence
When **AUDIO** key is pressed, deactivates audible alerts for two minutes. Indicated by illuminated  (2 Min LED).
- Audio Off
Feature user selectable. If enabled, press and hold **AUDIO** key for 3 seconds, and audible alarms will not activate. Indicated by flashing  (OFF LED).
- Battery Backed Trend Memory
Trend memory print of any 30 minutes, 2 hours, 8 hours 12 hours or 24 hours when used with the Seiko DPU-414, Seiko DPU-411, Hewlett-Packard ThinkJet, or Model 315 Printer.
- Analog (Recorder) Output Module—Optional
Provides analog output for strip chart applications at the following levels;
Oxygen Saturation value, 10mV/% (100% = 1 V)
Pulse Rate value, 4mV/bpm (250 bpm = 1 V)
Plethysmograph pulse waveform, 0-1V max (AGC)

- Serial (RS232) Data Output
 - Provides RS232 data interface compatible with;
 - Seiko DPU-414 Thermal Printer
 - Seiko DPU-411 Thermal Printer
 - Hewlett-Packard ThinkJet Printer
 - Novamatrix Model 315 Printer
 - Novamatrix Model 1260 Capnograph
 - Novamatrix Model 1010 Telemetry Central Station
 - SARACAP® monitor
 - RS232 computer interface
 - NOVACARD Memory Module
 - NOVACOM1 Interface
- Internal Real Time Clock
- Alert Bar

17











Accessories

Model 520A

17.1

Catalog No. Description

SuperBright™ SpO₂ SENSORS

8793-00	OxySnap™ Y-Sensor™ (use with OxySnap™ Extension Cable), 90 day warranty	
8744-00	OxySnap™ Finger Sensor (use with OxySnap™ Extension Cable), 1 yr. warranty	
8853-00	OxySnap™ Extension Cable (use with OxySnap™ sensors), 8 feet, 1 yr. warranty	
8898-00	OxySnap™ Extension Cable (use with OxySnap™ sensors), 12 feet, 1 yr. warranty	
8776-00	SuperBright™ Finger Sensor (10 ft sensor cable), 1 yr. warranty	
8791-00	SuperBright™ Y-Sensor™ (10 ft sensor cable), 90 day warranty	
9768-00	Finger Sensor (use with DB-9 Sensors), 6 feet, 1 yr. warranty	
9168-00	Finger Sensor (use with DB-9 Sensors), 3 feet, 1 yr. warranty	
9769-00	Y-Sensor™ (use with DB-9 Sensors), 6 feet, 90 day warranty	
9169-00	Y-Sensor™ (use with DB-9 Sensors), 3 feet, 90 day warranty	

SENSOR MANAGEMENT PLANS

Select a Finger Sensor or Y-Sensor™ Management Plan for each SuperBright™ Pulse Oximeter. The plan you select determines the length of coverage—36 or 60 months.








How the Plans Work: Included in each Plan are TWO sensors—one for immediate use, the other one for back-up. If a sensor becomes inoperative, place the backup sensor into use and return the inoperative sensor in the convenient prepaid mailer. A replacement sensor will be shipped within two business days of receipt of the inoperative sensor. This simple return/replacement method will be used for the entire warranty period, thereby, guaranteeing your costs and virtually eliminating sensor tracking hassles.

Warranty: For each Pulse Oximeter a plan is purchased for, the warranty on the monitor is also extended to the length of the Plan (a pre-contract inspection may be required). Replacement sensors provided under terms of the Plan shall carry the remaining Plan warranty—replacements do not extend the warranty.






8791-36 **Y-36 Plan** The Plan length is 36 months. Includes 9 boxes of any Y-Strip Taping Systems






Catalog No.	Description
8791-60	Y-60 Plan The Plan length is 60 months. Includes 15 boxes of any Y-Strip Taping Systems
8776-36	Finger-36 Plan The Plan length is 36 months
8776-60	Finger-60 Plan The Plan length is 60 months

Y-SENSOR™ APPLICATORS






8828-00	20mm Wrap Style Taping System (100 per box) Use on neonatal foot and hand, or on pediatric toe or finger, color coded blue	
8829-00	25mm Wrap Style Taping System (100 per box) Use on neonatal foot and hand, color coded green	
8831-00	20mm Finger Style Taping System (100 per box) Use on pediatric finger or on small adult finger, color coded blue	
8832-00	25mm Finger Style Taping System (100 per box) Use on adult finger, color coded green	
6929-00	Adhesive Foam Wraps, Large (25 per box)	
6968-00	Adhesive Foam Wraps, Small (25 per box)	
8836-00	Non-Adhesive Foam Wraps, Large (25 per box)	
8943-00	Non-Adhesive Foam Wraps, Small (25 per box)	
6131-50	Ear Clips (5 per box)	
6131-25	Ear Clips (25 per box)	

OUTPUT OPTIONS—ANALOG, DIGITAL & PRINTER

5963-00	Analog Output Module , for Model 2001	
9622-01	Analog Output Module , for OXYPLETH®, includes RS232 pass-thru	
6045-00	Cable , for 5963-00 Analog Output Module (open ended, 6 ft)	
5333-00	Cable , for 9622-00 Analog Output Module (open ended, 6 ft)	
2736-00	Analog/RS232 Output Option Kit for Model 515B (if ordered at time of purchase)	
6659-00	Cable , for 2736-00 Analog Output Option Kit (open ended, 6 ft)	

Catalog No.	Description	
5334-00	Cable, Serial Output to Personal Computer (PC with 25-pin connector)	
5335-00	Cable, Serial Output to Personal Computer (PC with 9-pin connector)	
5726-00	Cable, Analog Output Module (9622-01) to Aequitron 9500 Series Monitor	
9140-00	Thermal Printer, Seiko DPU-414 (with battery pack)	
400052	AC Adapter, for Seiko DPU-414 Printer, 120 VAC	
300017	Printer Paper, for Seiko DPU-411/414 Printer (2 rolls per box)	
9085-00	Printer Cable, for Seiko DPU-414 Printer, 9-to-25 pin (Models 515B/520A/ 860/1265/7100/2001)	

PULSE OXIMETRY ACCESSORIES

420034	Model 515B/C Inservice Video, VHS video tape (NTSC format)	
420016	OXYPLETH® Model 520A Inservice Video, VHS video tape (NTSC format)	
7106-10	Transport Pouch, for Models 515/520A/860/1265/7100/2001	
7104-10	Side Accessory Pouch, (included with monitor)	
600026	Power Cord, (included with monitor)	
4941-00	Saturation Sensor Extension Cable (4 feet)	
4942-00	Saturation Sensor Extension Cable (6 feet)	
4943-00	Saturation Sensor Extension Cable (10 feet)	
5266-00	Saturation Sensor Extension Cable (25 feet)	
6147-00	Saturation Sensor Extension Cable (50 feet)	
6455-00	Single Patient Use Pediatric/Adult Sensor (10 per box)	
6455-25	Single Patient Use Pediatric/Adult Sensor (25 per box)	
6480-00	Single Patient Use Neonatal/Pediatric Sensor (10 per box)	
6480-25	Single Patient Use Neonatal/Pediatric Sensor (25 per box)	
8933-00	Sensor Extension Cable, for use with Single Patient Use Sensors (DB-9 connectors)	
8936-00	Sensor Jumper Cable, for use with Single Patient Use Sensors (DB9 to OxySnap™ connector)	

EXTENDED WARRANTY

	Normal warranty: Monitor—1 year
9400-81	Model 2001 warranty extended an additional 1 year(s) at time of purchase
5693-81	OXYPLETH® warranty extended an additional 1 year(s) at time of purchase
6500-81	Model 515B warranty extended an additional 1 year(s) at time of purchase
6550-81	Model 515C warranty extended an additional 1 year(s) at time of purchase

BIOMEDICAL EQUIPMENT & TRAINING

Catalog No.	Description
	Service Test Kits include items and materials qualified service personnel may require to determine the functional integrity and/or accuracy of the system.
5777-00	Test Kit , Pulse Oximeter
9999-96	“focus” Technical Training Seminar (1 day course) For class schedules call: 1-800-243-3444 Ext. 2565

MOUNTING SYSTEMS

ROLLSTAND

140098	Rollstand , includes base, casters, tilt swivel head, post, poll, handle, utility basket. Requires either 140095 or 140101 Mounting Plate Kit.
--------	---

ROLLSTAND AND MOUNTING ACCESSORIES

140095	Adapter Plate Kit , for 515B/C, 520A, 2001, 1265/7100, 860, and <i>CO₂SMO Plus!</i> [®]
140099	Laptop Tray with Adapter Plate
140030	Wall Mount
140031	Wall Mount (less Wall Channel)
140032	Pivot Block Mount
140036	Countertop Mount (5 inch Base)
140100	Swivel C-Clamp Hanger
140067	Clamping Block
140068	Counterweight , for Rollstand
140070	Pole Mount Bracket , (for 2" pole), for Models 515B/C, 520A, 2001, 1265/7100, 860, and <i>CO₂SMO Plus!</i> [®]
140093	Mounting Plate Upgrade Kit , contains Mounting Plate Only
140094	Pole Mount , (for 3/4" - 2" dia. pole), for Models 515B/C, 520A, 2001, 1265/7100, 860, and <i>CO₂SMO Plus!</i> [®]
140097	6" Utility Basket , for Rollstand

[This page intentionally blank]


18

Parts Lists

Family Tree

18.1

The Model 520A diagram below shows the individual assemblies of the finished product—the 5693-00 Final Assembly.



5693-01 MAIN ASSEMBLY

18.2

Seq	Item Nbr	Description	Qty
013	140002	CARD GUIDE, 2.5" L	2
005	5713-01	SPEAKER ASSY	1
007	5719-01	FRONT PANEL ASSY	1
014	161067	TAPE, CL CELL	0

Seq	Item Nbr	Description	Qty
010	5766-01	BOTTOM COVER ASSY	1
008	5728-01	CABLE ASSY, SPO2 IN	1
001	2471-01	MAIN BOARD ASSY,	1
012	9621-16	LENS, ALERT	1
017	286220	SCREW, 6-32 X 3/16	0
003	6538-01	REAR PANEL ASSY	1
002	5664-10	CHASSIS, MODEL 520A	1
016	286205	6-32 x 3/4 S.B.H.	0
006	5714-01	BATTERY HARNESS ASSY	1
015	284200	#4-40 x 1/4 SLOTTED	0
018	400024	BATTERY, 12VDC, 2.3A	1
009	5765-01	TOP COVER ASSY	1
001	2775-01	MAIN BOARD ASSY	1

6538-01 REAR PANEL ASSY

18.3

Seq	Item Nbr	Description	Qty
007	6538-17	REAR PANEL SUBASSY	1
018	285001	LOCK WASHER, NO. 6	0
019	285005	FLAT WASHER, NO. 8	0
003	5812-10	SPACER, SUPPORT	2
002	4109-10	SUPPORT BRACKET, PC	1
024	161092	ADHESIVE, 262	0
006	9090-01	POWER CABLE ASSY	1
001	2726-01	POWER SUPPLY BD ASSY	1
014	280188	STANDOFF, 3/8 DIA	0
004	5820-01	GROUND WIRE ASSY	1
012	210149	POWER ENTRY MODULE,	1
023	285013	NYL WASH, #4 NATURAL	0
017	285000	LOCK WASHER, NO. 4	0
022	608033	WIRE CLIP & BUSHING	1
013	216059	CONNECTOR, PLUG	1
011	161008	ADHESIVE, 242	0
021	515023	FUSE, 1/2A, 250V	2
020	286219	SCREW, 6-32 X 1 1/4	0
008	6549-01	GROUND WIRE ASSY	1
016	281501	NUT, HEX, NO. 6-32	0
015	281500	NUT, HEX, NO. 4-40	0

5673-01 REAR PANEL ASSY

18.4

Seq	Item Nbr	Description	Qty
016	285001	LOCK WASHER, NO. 6	0
007	5812-10	SPACER, SUPPORT	2
002	4109-10	SUPPORT BRACKET	1
019	5819-01	GROUND WIRE ASSY	1
022	285046	WASHER, FLAT, NO. 6	0

Seq	Item Nbr	Description	Qty
005	5673-32	REAR FACEPLATE	1
006	5743-10	POWER ENTRY MODULE	1
013	280188	STANDOFF, 3/8 DIA	0
020	5820-01	GROUND WIRE ASSY	1
003	5229-11	WIRE CLIP MOD	1
017	285013	NYL WASH, #4 NATURAL	0
015	285000	LOCK WASHER, NO. 4	0
011	216059	CONNECTOR, PLUG	1
023	161008	ADHESIVE, 242	0
008	9308-11	BUSHING MOD, LINECORD	2
021	515023	FUSE, 1/2A, 250V	2
018	286219	SCREW, 6-32 X 1 1/4	0
004	5673-10	REAR PANEL	1
012	281501	NUT, HEX, NO. 6-32	0
001	2472-01	POWER SUPPLY BOARD	1
014	281500	NUT, HEX, NO. 4-40	0

2726-01 POWER SUPPLY BD ASSY**18.5**

Seq	Item Nbr	Description	Qty
000	211505	CONNECTOR, 5 PIN	1
000	474145	RESISTOR, 215 OHM	
000	481530	DIODE, KBU4G, 4 AMP	
000	474214	RESISTOR, 332K OHM	1
000	474137	RESISTOR, 1M OHM	1
000	481549	DIODE, MBRS14OT3	1
000	152096	CAPACITOR, 220UF	1
000	484529	VOLTAGE REGULATOR	1
000	180014	INDUCTOR, 2SUH	1
000	474165	RESISTOR, 10K OHM	1
000	152081	CAPACITOR, 6BOUF	1
000	474162	RESISTOR, 61.9K OHM	1
000	5918-10	XFMR, MAINS, PWR SPL	1
000	470026	RESISTOR, 150 OHM	1
000	474215	RESISTOR, 2.21K OHM	1
000	486805	IC, LM393M, DUAL VOL	1
000	2725-02	RAE, POWER SUPPLY BD	1
000	515085	FUSE W FUSEHOLDER	1
000	474141	RESISTOR, 249K OHM	1
000	2726-04	TEST PROCEDURE, PWR	0
000	600034	RIBBON CABLE ASSY	1
000	474181	RESISTOR, 4.3 OHM	1
000	474166	RESISTOR, 100K OHM	4
000	474216	RESISTOR, 4.99K OHM	1
000	485543	TRANSISTOR, MOSFET	1
000	211213	CONNECTOR, 2 PIN	1
000	481552	DIODE, MBRS34OT3	3
000	154072	CAPACITOR, .1UF	4

Seq	Item Nbr	Description	Qty
000	474151	RESISTOR, 37.4K OHM	1
000	484559	VOLT RGLR, LT1O76CT	1
000	515083	FUSE W FUSEHOLDER	1
000	2726-17	POWER SUPPLY BOARD	1
000	154079	CAP, 1OUF, 2EV, 20%,	1
000	474218	RESISTOR, 47.5K OHM	1
000	280114	SNAP RIVET, .118	0
000	152029	CAPACITOR, 220UF 20%	1
000	485532	TRANSISTOR, 2N7002T1	1
000	474211	RESISTOR, 49.9K OHM	1
000	481561	DIODE, ULTRAFASST, 3A	1
000	474341	RESISTOR, 4.3 OHM	1
000	481563	DIODE, SCHOTTKY, 5A	2

2472- 01 POWER SUPPLY BOARD**18.6**

Seq	Item Nbr	Description	Qty
024	472065	RESISTOR, 200K OHM	1
016	211505	CONNECTOR, 5 PIN	1
036	484539	VOLTAGE REGULATOR	1
001	2472-02	FAB, POWER SUPPLY BD	1
032	481530	DIODE, KBU4G, 4 AMP	1
040	515046	FUSE, 2A, 250V	1
020	472022	RESISTOR, 4.99K OHM	1
044	161049	ADHESIVE, HOT MELT	0
028	472185	RESISTOR, 61.9K OHM	1
035	484529	VOLTAGE REGULATOR	1
014	180014	INDUCTOR, 25UH, .25	1
027	472140	RESISTOR, 332K OHM	1
012	152081	CAPACITOR, 6800UF	1
037	485517	TRANSISTOR, VN10KM	1
041	515035	FUSE, 1/2A, 250V	1
034	481542	DIODE, UF5400	2
033	481541	DIODE, 100V, 1A	1
004	5918-10	XFMR, MAINS, PWR SPL	1
017	470026	RESISTOR, 150 OHM	1
023	472058	RESISTOR, 100K OHM	3
018	472005	RESISTOR, 249 OHM	1
043	600034	RIBBON CABLE ASSY	1
030	474181	RESISTOR, 4.3 OHM	1
039	486716	LM393N DUAL COMPARATOR	1
022	472049	RESISTOR, 49.9K OHM	1
019	472016	RESISTOR, 2.21K OHM	1
038	485527	MOSFET, IRF9523	1
031	481022	DIODE, ZENER, BZY88C	1
026	472146	RESISTOR, 47.5K OHM	1
042	515503	FUSE HOLDER, FOR 5MM	4
021	472030	RESISTOR, 10K OHM	1

Seq	Item Nbr	Description	Qty
015	211213	CONNECTOR, 2 PIN	1
013	154016	CAPACITOR, .1UF	5
011	152029	CAPACITOR, 220UF	1
025	472109	RESISTOR, 37.4K OHM	1
029	472270	RESISTOR, 1.2M OHM	1

5719-01 FRONT PANEL ASSY**18.7**

Seq	Item Nbr	Description	Qty
006	280033	SPACER *4 X 1/8	0
003	5720-01	DISPLAY ASSY	1
001	2473-01	ALERT BOARD ASSY	1
010	285000	LOCK WASHER, NO. 4	0
002	5676-27	MEMBRANE KEYANEL	1
008	284200	*4-40 X 1/4 SLOTTED	0
007	280187	STANDOFF, .187 DIA	0
004	6278-10	SHIELD, VFD	1
009	284204	*4-40 X 1/2 SLOTTED	0

2473-01 ALERT BOARD ASSY**18.8**

Seq	Item Nbr	Description	Qty
006	482565	LED, CLEAR, ROUND	4
005	470024	RESISTOR, 100 OHM	2
007	600038	RIBBON CABLE ASSY	1
001	2473-02	FAB, ALERT BOARD	1

5720- 01 DISPLAY ASSY**18.9**

Seq	Item Nbr	Description	Qty
001	482573	LCD MODULE, DOT MATRIX	1
002	600033	RIBBON CABLE ASSY	1

5713-01 SPEAKER ASSY**18.10**

Seq	Item Nbr	Description	Qty
001	130010	SPEAKER, 8 OHM 15%,	1
002	211225	CONNECTOR, 2 PIN	1
005	5855-10	FOAM LINER	1
004	608003	TUBING, HEAT SHRINK	0
003	605060	WIRE, 26AWG, PVC INS	0

5714-01 BATTERY HARNESS ASSY**18.11**

Seq	Item Nbr	Description	Qty
004	605163	WIRE, 22AWG, BLACK	0

Seq	Item Nbr	Description	Qty
001	211215	CONNECTOR, 2 PIN	1
002	216068	TERMINAL LUG, RECEPT	0
005	608001	CABLE TIE	0
003	605162	WIRE, 22AWG, RED	0

5728-01 CABLE ASSY

18.12

Seq	Item Nbr	Description	Qty
004	211704	CONNECTOR, RCPT	1
011	605063	WIRE, 26AWG, PVC INS	0
015	608005	TUBING, HEAT SHRINK	0
006	210098	CONTACT, FEMALE	0
008	605059	WIRE, 26AWG, PVC INS	0
007	280108	SNAP RIVET	0
001	5606-16	DRESS BEZEL, SAO2 IN	1
014	608001	CABLE TIE	0
010	605061	WIRE, 2EAWG, PVC INS	0
003	180009	FERRITE SLEEVE	1
013	605065	WIRE, 26AWG, PVC INS	0
012	605064	WIRE, 26AWG, PVC INS	0
002	161060	EPOXY, 2 PART (5/7)	0
005	211630	CONNECTOR, 6 PIN	1
009	605060	WIRE, 26AWG, PVC INS	0
016	608011	TUBING, HEAT SHRINK	0
002	161110	EPOXY, 2 PART (1/1)	0
017	9779-33	ASSY PROC, EPOXY MIX	0

5765-01 TOP COVER ASSY

18.13

Seq	Item Nbr	Description	Qty
005	5862-10	SHIELD, MYLAR	1
001	5765-13	TOP COVER	1
002	5827-32	WARNING LABEL	1
004	5823-32	LABEL, INSTRUCTIONS	1
003	5828-32	LABEL, ISOLATION	1

5766-01 BOTTOM COVER ASSY

18.14

Seq	Item Nbr	Description	Qty
016	161007	ADHESIVE	0
010	284261	SCREW, 4-40 X 5/8 L	0
009	4727-10	KICKSTAND, BEDRAIL	1
005	5826-10	FOOT PAD	4
003	5766-13	BOTTOM COVER	1
015	284264	SCREW, 4-40 X 1/4 L	0
011	161064	TAPE, 3/4 X 60 YDS	0
012	5405-10	SHIELD, BATTERY	1

Seq	Item Nbr	Description	Qty
007	5760-16	LEFT FOOT, WHITE	2
008	5761-16	RIGHT FOOT, WHITE	2
013	5849-10	BRACKET, SUPPORT	1
002	5409-32	LABEL, WARNING, POLE	1
006	315052	LABEL, EARTHING SYMB	13

2755-01 TONE GEN REPL BOARD**18.15**

Seq	Item Nbr	Description	Qty
022	486358	IC, NC7SZ32MS, 2-IN	1
018	474241	RESISTOR, 150K OHM	1
023	486825	IC, LMC7101BIM5X	1
020	474300	RESISTOR, 402K OHM	1
014	474228	RESISTOR, 1.21K OHM	1
001	2755-02	FAB, TONE GENERATOR	1
017	474240	RESISTOR, 100K OHM	3
012	474197	RESISTOR, 49.9K OHM	1
013	474224	RESISTOR, 100 OHM	1
011	211814	CONNECTOR, 8 PIN	2
010	154104	CAPACITOR, .01UF	4
002	2755-03	SCHEM, TONE GEN REPL	0
003	9074-07	IC, TONE GENERATOR	1
016	474235	RESISTOR, 24.9K OHM	2
019	474270	RESISTOR, 200K OHM	1
011	211815	CONN, 8 PIN, HDR	2

2775-01 MAIN BOARD ASSY**18.16**

Seq	Item Nbr	Description	Qty
000	472121	RESISTOR, 3.74K OHM	1
000	483017	TRANSISTOR, BC214C	4
000	474032	RESISTOR, 10 OHM	1
000	472041	RESISTOR, 20K OHM	2
000	472063	RESISTOR, 162K OHM	1
000	472200	RESISTOR, 5.6 OHM	1
000	487064	IC, CNY17-III, PHOTO	4
000	154058	CAPACITOR, .022UF	1
000	472011	RESISTOR, 1K OHM	5
000	474089	RESISTOR PACK, 100K	2
000	474131	RESISTOR PACK, 10K	1
000	211306	CONNECTOR, 3 PIN	1
000	180033	EMI FILTER, 22PF	1
000	472255	RESISTOR, 1 OHM	2
000	154065	CAP, .1UF, 63V	5
000	152040	CAPACITOR, 2.2UF	3
000	472007	RESISTOR, 475 OHM	6
000	154081	CAPACITOR, 100PF	5

Section 18 *Parts Lists*

Seq	Item Nbr	Description	Qty
000	472021	RESISTOR, 4.75K OHM	1
000	210144	JUMPER, 2 POSITION	1
000	153003	CAPACITOR, .01UF	1
000	472193	RESISTOR, 23.7K OHM	1
000	481541	DIODE, bOy, 1A	4
000	475034	POTENTIOMETER, 2K OHM	1
000	483002	2N3906 PNP SILICON	1
000	472058	RESISTOR, 100K OHM	7
000	154060	CAPACITOR, .22UF,	1
000	481501	DIODE, 1N4148, SIGNAL	9
000	472274	RESISTOR, 26.7 OHM	1
000	472105	RESISTOR, 3.92K OHM	2
000	470109	RESISTOR, 470K OHM	1
000	481534	DIODE, BAT82, SCHOTT	2
000	2471-28	MAIN BD KIT	1
000	472227	RESISTOR, 47 OHM	2
009	285000	LOCK WASHER, NO. 4	0
000	474134	RESISTOR PACK, 41K	1
000	486299	IC, CS5SO3JP, 20-BIT	2
000	472170	RESISTOR, 22.1K OHM	1
000	474132	RESISTOR PACK, 41K	2
000	472276	RESISTOR, 6.49K OHM	2
000	470016	RESISTOR, 22 OHM	1
000	472268	RESISTOR, 8.66K OHM	1
000	5833-01	GROUND WIRE ASSY	1
000	472003	RESISTOR, 100 OHM	2
000	472146	RESISTOR, 47.5K OHM	1
000	152072	CAPACITOR, 10UF	1
000	472246	RESISTOR, 66.5K OHM	2
000	472198	RESISTOR, 3.3 OHM	1
000	474098	RESISTOR PACK, 100K	4
000	153027	CAPACITOR, 220PF	1
000	472030	RESISTOR, 10K OHM	7
000	474133	RESISTOR PACK, 10 OH	1
000	474086	RESISTOR PACK, 10K	1
000	472034	RESISTOR, 12.1K OHM	1
000	153063	CAPACITOR, 220PF	2
000	152075	CAPACITOR, 47UF, 25V	2
000	474138	RESISTOR, 100 OHM,	4
000	154016	CAPACITOR, .1UF	59
000	216029	TEST POINT, SPRING	2
000	152045	CAPACITOR, 10UF 16V	6
000	153052	CAPACITOR, .022UF	1
000	152073	CAPACITOR, 100UF	4
007	280023	SPACER HEX #4-40	0
000	153021	CAPACITOR, .47UF	1
000	472195	RESISTOR, 41.2K OHM	4
000	180011	FERRITE BEAD, 22 AWG	1

Seq	Item Nbr	Description	Qty
000	153051	CAPACITOR, .22UF,	1
000	474135	RES PACK, 100 OHM	1
000	481031	DIODE, ZENER, BZX79-	1
001	2755-01	TONE GEN REPL BOARD	1
000	472037	RESISTOR, 13.7K OHM	2
008	284217	SCREW, 4-40 X 3/16	0
000	180067	FERRITE BEAD, 600 uH	15
000	215080	SOCKET, 32 PIN, PLCC	1
000	481050	ESD SUPPRESSOR	2
002	2775-02	FAB, MAIN BD	1
000	212136	CONN, 10 PIN, HDR	1
000	5715A-07	PROGRAM, SPO2, 520A	1
000	211237	CONNECTOR, 2 PIN	1
000	486845	IC, 2.SV OR 3.OV V R	1
006	161039	FOAM TAPE, 1/32THK	0
010	5712-10	TRANSFORMER, MAIN BD	1
010	5711-10	TRANSFORMER, MAIN BD	0

2471-28 MAIN BD KIT (part of 2775-01)**18.17**

Seq	Item Nbr	Description	Qty
001	180012	FERRITE BEAD	1
000	484535	VOLTAGE REGULATOR, L	1
000	486295	IC, MM74HC40S1N	1
000	484534	VOLTAGE REGULATOR	1
000	486625	MC74HC32N I.C. QUAD	3
000	211414	CONNECTOR, 4 PIN	1
000	152085	CAPACITOR, 1500UF	2
000	515511	FUSEHOLDER, VERT	2
000	484533	VOLTAGE RGLTR	1
000	484523	VOLTAGE REGULATOR	1
000	484529	VOLTAGE REGULATOR	1
000	487061	IC, TBA820M, MONOLIT	1
000	486300	IC, TLC549IP, 8-BIT	1
000	487069	IC, LT1019CN8-2.5	1
000	486712	TLO74CN QUAD F FET	2
000	152086	CAPACITOR, 1000UF	2
000	485528	MOSFET, N-CH	1
000	486651	MM74HC138 I.C.	1
000	471400	RESISTOR, IOOM OHM	1
000	153006	CAPACITOR, 47PF, 63V	4
000	486298	IC, MM74HC402ON	1
000	486680	IC, MM74HC573AN	4
000	9715-07	PROGRAM, PEEL ASSY	1
000	9716-07	PROGRAM, PEEL ASSY	1
000	154076	CAPACITOR, .01UF	2
000	481542	DIODE, UF5400	1
000	215060	SOCKET, IC, 32 PLN	1

Section 18 *Parts Lists*

Seq	Item Nbr	Description	Qty
000	230006	CRYSTAL,3.2768MHZ	1
000	152084	CAPACITOR, 470UF	2
000	230018	CRYSTAL, 12.288MHZ	1
000	486268	IC, DG444DJ, SPST	1
000	486305	IC, MSM6242BRS, CMOS	1
002	153013	CAPACITOR, 33PF	3
02A	153013	CAPACITOR, 33PF	1
000	210097	CONNECTOR, 16 PIN	1
000	230016	CRYSTAL, 32.768K HZ	1
000	486600	IC, CD4013B	1
000	484531	VOLTAGE REGULATOR	1
000	486285	IC, HD64180R1P6	1
000	515071	FUSE, 1/10A, 125V	1
000	180010	INDUCTOR, 18UH	1
003	153012	CAPACITOR, 22PF	4
03A	153012	CAPACITOR, 22PF	6
000	212501	CONNECTOR, 20 PIN	1
000	486675	MM74HC14N, IC,CMOS	2
000	486685	IC, DUAL 8-BIT	1
000	486256	IC, K6TO8O8C1D-DB7O	1
000	485529	TRANSISTOR, VNO61OL	7
000	154057	CAPACITOR, 470pF	1
000	485527	MOSFET, 1RF9523	1
000	211412	CONNECTOR, 4 PIN	1
000	484014	NPN TRANSISTOR ZTX1O	1
000	486276	IC, LT1O81CN	1
000	180004	CHOKE, 100UH, 10%	1
000	215031	SOCKET, 20 PIN, DIP	2
000	211213	CONNECTOR, 2 PIN	2
000	211629	CONNECTOR, 6 PIN	1
000	470111	RESISTOR, .33 OHM	3
000	153046	CAPACITOR, 4700PF	1
000	486606	CD4O93BE	1
000	152066	CAPACITOR, 220uF	1
000	485520	TRANSISTOR, BS2SO	2
000	212529	CONNECTOR, 20 PIN	1
000	484515	LM317LZ REGULATOR	1
000	215055	SOCKET, IC	1
000	515072	FUSE, 1A, 125V	1
000	487065	IC, UC3843N, CURRENT	1
000	400035	INJERTER, DC TO AC	1
000	210051	CONNECTOR, 25 PIN	1
000	180011	FERRITE BEAD, 22 AWG	1
004	487053	IC, AD712JN, DUAL BI	2
004	486717	LF353N	2

2471-01 MAIN BOARD ASSY**18.18**

Seq	Item Nbr	Description	Qty
067	472121	RESISTOR, 3.74K OHM	1
153	483017	TRANSISTOR, BC214C	4
082	474032	RESISTOR, 10 OHM	1
062	472041	RESISTOR, 20K OHM	2
064	472063	RESISTOR, 162K OHM	1
073	472200	RESISTOR, 5.6 OHM	1
130	487064	IC, CNY17-III, PHOTO	4
151	154058	CAPACITOR, .022UF	1
057	472011	RESISTOR, 1K OHM	5
083	474089	RESISTOR PACK, IOOK	2
009	6344-07	PROGRAM, EPROM ASSY	0
085	474131	RESISTOR PACK, 10K	1
149	472255	RESISTOR, 1 OHM	2
032	154065	CAP, .1UF, 63V	5
150	486255	IC, SN76496AN	1
011	152040	CAPACITOR, 2.2UF	3
056	472007	RESISTOR, 475 OHM	6
002	2471-03	SCHEMATIC, MAIN BD	0
058	472021	RESISTOR, 4.75K OHM	1
004	5715-07	PROGRAM, SPO2, 520A	1
001	2471-02	FAB, MAIN BOARD	1
020	153003	CAPACITOR, .01UF	1
070	472193	RESISTOR, 23.7K OHM	1
094	481541	DIODE, bOy, 1A	4
090	475034	POTENTIOMETER, 2K	1
096	483002	2N3906 PNP SILICON	1
063	472058	RESISTOR, IOOK OHM	7
155	153062	CAPACITOR, 150PF	3
152	154060	CAPACITOR, .22UF	1
092	481501	DIODE, 1N4148, SIGNAL	10
080	472274	RESISTOR, 26.7 OHM	1
066	472105	RESISTOR, 3.92K OHM	2
079	472272	RESISTOR, 59K OHM	1
093	481534	DIODE, BAT82, SCHOTT	2
008	2471-28	MAIN BD KIT	1
075	472227	RESISTOR, 47 OHM	2
088	474134	RESISTOR PACK, 41K	1
003	5712-10	TRANSFORMER, MAIN BD	1
117	486299	IC, CSS5O3JP, 20-BIT	2
069	472170	RESISTOR, 22.1K OHM	1
086	474132	RESISTOR PACK, 41K	2
081	472276	RESISTOR, 6.49K OHM	2
077	470016	RESISTOR, 22 OHM	1
078	472268	RESISTOR, 8.66K OHM	1
138	5833-01	GROUND WIRE ASSY	1
055	472003	RESISTOR, 100 OHM	4

Section 18 *Parts Lists*

Seq	Item Nbr	Description	Qty
068	472146	RESISTOR, 47.5K OHM	1
014	152072	CAPACITOR, 10UF	1
076	472246	RESISTOR, 66.5K OHM	2
072	472198	RESISTOR, 3.3 OHM	1
010	6266-57	PROGRAM, DISK ASSY	0
084	474098	RESISTOR PACK, 100K	4
148	153027	CAPACITOR, 220PF	1
059	472030	RESISTOR, 10K OHM	6
087	474133	RESISTOR PACK, 10 OHM	1
145	474086	RESISTOR PACK, 10K	1
003	5711-10	TRANSFORMER, MAIN BD	0
060	472034	RESISTOR, 12.1K OHM	1
025	153038	CAPACITOR, .001UF	1
016	152075	CAPACITOR, 47UF, 25V	2
028	154016	CAPACITOR, .1UF	59
142	216029	TEST POINT	2
012	152045	CAPACITOR, 10UF 16v	6
144	153052	CAPACITOR, .022UF	1
015	152073	CAPACITOR, 100UF	4
024	153021	CAPACITOR, .47UF	1
071	472195	RESISTOR, 41.2K OHM	4
027	153051	CAPACITOR, .22UF	1
089	474135	RES PACK, 100 OHM	1
091	481031	DIODE, ZENER	1
154	2755-01	TONE GEN REPL BOARD	1
061	472037	RESISTOR, 13.7K OHM	2

19


Schematic and Assembly Drawings

Drawing	Description
5693-09	Overall Wiring Diagram
2726-01	PSU Board Assy
2726-03	PSU Board Schematic
2472-01	PSU Board Assy
2472-03	PSU Board Schematic
5693-01	Oximeter Unit Assy
2471-01	Main Board Assy
2471-03	Main Board Schematic
2775-01	Main Board Assy
2775-03	Main Board Schematic
2755-01	Tone Generator Replacement Bd
2755-03	Tone Generator Replacement Bd
5728-01	Cable Assy, SpO2 Input
2473-03	Alert Lamp Board
5676-27	Keypanel Schematic
5693-00	Main Assy, Model 520A

	Optional Analog Module
9622-01	Analog/RS232 Module
2289-01	Analog/RS232 Bd. Assy
2289-03	Analaog/RS232 Bd. Schematic

Table 5. Schematic and Assembly Drawings


This document contains information which is proprietary and the property of Novamatrix Medical Systems Inc. and may not be reproduced, used or disclosed in whole or in part without the written consent of Novamatrix Medical Systems Inc.



EMI/RFI SHIELDING

DO NOT SCALE UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES (mm) BREAK ALL SHARP EDGES				TITLE	
04	N557	27May98	DRG	~	OVERALL WIRING DIAGRAM, OXYPLETH - 520A
03	N547	10Feb98	MFG	~	NOVAMETRIX MEDICAL SYSTEMS INC. WALLINGFORD, CT U.S.A. 06492
02	N377	19Apr96	ENGR	~	
01	N375	16Nov95	ENGR	~	
REV	R. NO.	DATE	FINISH	~	USED ON: D5693-00
				SIZE	D
				DRAWING NO.	5693
				CODE	09
				REV	04
				SCALE:	NONE
				SHEET	1 OF 1


This document contains information which is proprietary and the property of Novamatrix Medical Systems Inc and may not be reproduced, used or disclosed in whole or in part without the written consent of Novamatrix Medical Systems Inc.



DETAIL A

REV	R NO.	DATE	REV	R NO.	DATE	FINISH	TITLE	NOVAMATRIX MEDICAL SYSTEMS INC. WALLINGFORD, CT U.S.A. 06492											
13	N727	17Jan00	07	N375	13Nov95	DO NOT SCALE UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES (mm)	MAIN ASSEMBLY, MODEL 520A	<table border="1"> <tr> <td>SIZE</td> <td>DRAWING NO.</td> <td>CODE</td> <td>REV</td> </tr> <tr> <td>D</td> <td>5693</td> <td>01</td> <td>13</td> </tr> </table>	SIZE	DRAWING NO.	CODE	REV	D	5693	01	13			
SIZE	DRAWING NO.	CODE	REV																
D	5693	01	13																
12	N557	26May98	06	N366	19May95	BREAK ALL SHARP EDGES TOLERANCES (mm) ± ~ FRAC ± ~ INCHES +.007 -.003 (mm) (±.18 -.08)	<table border="1"> <tr> <td>DRAWN</td> <td>CHK</td> <td>CHECKED</td> </tr> <tr> <td>MFL</td> <td>EH</td> <td>MRL</td> </tr> <tr> <td>WFO</td> <td>ENGR</td> <td>APPROVED</td> </tr> <tr> <td>AK</td> <td>22May92</td> <td></td> </tr> </table>	DRAWN	CHK	CHECKED	MFL	EH	MRL	WFO	ENGR	APPROVED	AK	22May92	
DRAWN	CHK	CHECKED																	
MFL	EH	MRL																	
WFO	ENGR	APPROVED																	
AK	22May92																		
11	N547	12Feb98	04	N325	15Jul93	MATERIAL	USED ON: D5693-00	SCALE: 1/1											
10	N427	10Jan97	03	N316	24Mar93														
09	N418	20Dec96	02	N304	15Oct92														
08	N377	10Jan96	01	N302	26May92														

This document contains information which is proprietary and the property of Novametric Medical Systems Inc. and may not be reproduced, used or disclosed in whole or in part without the written consent of Novametric Medical Systems Inc.




PN	DESCRIPTION	REF DES
		1
2725-02	FAB, POWER SUPPLY BOARD, 2ND GENERATION	~
D2726-03	SCHEMATIC, POWER SUPPLY BOARD	~
A2726-04	TEST PROCEDURE, POWER SUPPLY BOARD	~
☆D5918-10	TRANSFORMER, MAINS, POWER SUPPLY BOARD	T301
☆152029	CAPACITOR, 220uF 20%, 35V, ELECTROLYTIC	C4
☆152081	CAPACITOR, 6800uF, 35V, 20%, ELECTROLYTIC	C1
152096	CAPACITOR, 220uF, 35V, 20%, ELECTROLYTIC	C2
154072	CAPACITOR, .1uF, 50V, 10%, X7R, CERAMIC CHIP	C3, C5-C7
154079	CAPACITOR, 10uF, 25V, 10%, TANTALUM	CB
☆180014	INDUCTOR, 25uH, .25 SPACING	L1
☆211213	CONNECTOR, 2 PIN, POST HEADER	J302
☆211505	CONNECTOR, 5 PIN, SQUARE POST HEADER	J301
280114	SNAP RIVET, .118--.158 THK PANEL	FOR IC1
☆470026	RESISTOR, 150 OHM, 1/4W, 10%	R15
474137	RESISTOR, 1M OHM, 1/8W, 1%	R1
474141	RESISTOR, 249k OHM, 1/8W, 1%	R2
474145	RESISTOR, 215 OHM, 1/8W, 1%	R11
474151	RESISTOR, 37.4k OHM, 1/8W, 1%	R6
474162	RESISTOR, 61.9k OHM, 1/8W, 1%	R13
474165	RESISTOR, 10k OHM, 1/8W, 1%	R10
474166	RESISTOR, 100k OHM, 1/8W, 1%	R5, R7, R9, R17
474211	RESISTOR, 49.9k OHM, 1/8W, 1%	R4
474214	RESISTOR, 332k OHM, 1/8W, 1%	R14
474215	RESISTOR, 2.21k OHM, 1/8W, 1%	R16
474216	RESISTOR, 4.99k OHM, 1/8W, 1%	R12
474218	RESISTOR, 47.5k OHM, 1/8W, 1%	R8
☆474341	RESISTOR, 4.3 OHM, 3W, 5%	R3
☆481530	DIODE, 4 AMP, BRIDGE RECTIFIER	D1
481549	DIODE, RECTIFIER	D2
481561	DIODE, ULTRAFAST, 100V, 3A	D4
481563	DIODE, SCHOTTKY, 100V, 5A	D3, D5
☆484529	VOLTAGE REGULATOR, 5V, 100MA	IC2
☆484559	VOLTAGE REGULATOR, FLOW 30, 2A	IC1
485532	TRANSISTOR, N-CHANNEL, ENHAN MODE	Q2
485543	TRANSISTOR, MOSFET, N-CHANNEL & P-CHANNEL	Q1
486805	IC, DUAL VOLTAGE COMP, LOW POWER	IC3
515083	FUSE W FUSEHOLDER, 1A, 125V, V FAST-ACTING	F302
515085	FUSE WITH FUSEHOLDER, 2A, 125V, SLO-BLO	F301
600034	RIBBON CABLE ASSY, 4 PIN, 4 INCHES LONG	P102

NOTES:
 1. ☆ DENOTES THRU HOLE COMPONENTS.
 2. FOR TEST PROCEDURE SEE A2726-04.

REV	R NO.	DATE	FINISH	DO NOT SCALE UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES (mm) BREAK ALL SHARP EDGES TOLERANCES (mm) ± DEC ± .007 (mm) ± .003 (mm) (+.18 -.08) HOLES ± .007 (mm) (+.18 -.08) MATERIAL	TITLE	DRAWN	CHKD	CHECKED	BL	SIZE	DRAWING NO.	CODE	REV
06	N675	28Jun99			POWER SUPPLY BOARD ASSY, MODEL 515A SB, 520A, 1265 & 7100					D	2726	01	06
05	N671	1Jun99											
04	N616	3Nov98											
03	N508	29Aug97											
02	N433	29Jan97											
01	N377	26Dec95											
00	N372	29Sep95											
REV	R NO.	DATE	FINISH										




WALLINGFORD, CT U.S.A. 06492



- NOTES:
- UNLESS OTHERWISE SPECIFIED RESISTANCE VALUES ARE IN OHMS, 1/8 W, 1%.
 - IC1 IS THERMALLY COUPLED TO THE CIRCUIT BOARD.
 - CONNECT DGND AND AGND ONLY AT IC1-3
 - 2725 SPECIFICATIONS:
 - A. COMPONENTS NOT FITTED: D3, F301, IC2, R15, P102
 - B. F302: 2A, SLO-BLO
 - C. R6: 26.1k
 - 2726 SPECIFICATIONS:
 - A. F302: 1A, VERY FAST-ACTING
 - B. R6: 37.4k
 - LAST REF DES USED: C8, D5, E302, F302, IC3, J302, L1, Q3, R17

SCHEMATIC, POWER SUPPLY BOARD				Novametric Medical Systems, Inc. 1 Barnes Industrial Park Road Wallingford, CT USA 06492			
06	N675	28Jun99					
05	N671	21Jun99					
04	N616	13Nov98	Drawn: R. Daniels	Checked: MRL	30Mar95	3Oct95	
03	N508	2Sept97	Mfg E. Harvie	Approved: AJE	7Nov95	7Nov95	
02	N433	27Feb97	Eng: 7Nov95				
01	N377	26Dec95					
00	N372	3Oct95	File: 2725-260306				
Rev.	"R" No.	By	Date	Date: Monday, June 28, 1999			

This document contains information which is proprietary and the property of Novamatrix Medical Systems Inc. and may not be reproduced, used or disclosed in whole or in part without the written consent of Novamatrix Medical Systems Inc.




PN	DESCRIPTION	REF DES
5918-10	TRANSFORMER, MAINS	T301
152029	CAPACITOR, 220uF, 35V, ELECTROLYTIC	C4
152081	CAP, 6800uF, 35V, ±20%, AL, ELECTROLYTIC	C1
154016	CAPACITOR, .1uF, 50V, DIPPED CERAMIC	C2, C3, C5-C7
180014	INDUCTOR, 25uH	L1
211213	CONNECTOR, 2 PIN, POST HEADER	J302
211505	CONNECTOR, 5 PIN, POST HDR, FRICTION LOCK	J301
470026	RESISTOR, 150 OHM, 1/4W, ±10%	R15
472005	RESISTOR, 249 OHM, 1/4W, ±1%	R11
472016	RESISTOR, 2.21k OHM, 1/4W, ±1%	R16
472022	RESISTOR, 4.99k OHM, 1/4W, ±1%	R12
472030	RESISTOR, 10k OHM, 1/4W, ±1%	R10
472049	RESISTOR, 49.9k OHM, 1/4W, ±1%	R4
472058	RESISTOR, 100k OHM, 1/4W, ±1%	R5, R7, R9
472065	RESISTOR, 200k OHM, 1/4W, ±1%	R1
472109	RESISTOR, 37.4k OHM, 1/4W, ±1%	R8
472146	RESISTOR, 47.5k OHM, 1/4W, ±1%	R8
472140	RESISTOR, 332k OHM, 1/4W, ±1%	R14
472185	RESISTOR, 61.8k OHM, 1/4W, ±1%	R13
472270	RESISTOR, 1.2M OHM, 1/4W, ±1%, METAL FILM	R2
474181	RESISTOR, 4.3 OHM, 1W, ±5%	R3
481022	DIODE, ZENER, BZY88C12, 12V, 500mW	D2
481530	DIODE, KBU4D, BRIDGE RECTIFIER, 4 AMP	D1
481541	DIODE, UF4002, 100V, 1A, FAST RECOVERY	D5
481542	DIODE, UF5400, 50V, 3A, FAST RECOVERY	D3, D4
484529	VOLTAGE REGULATOR, LM78L05ACZ, 5V, 100mA	IC2
484539	VOLTAGE REGULATOR, LT1076CT, 2A, TO-220	IC1
485517	TRANSISTOR, VN10KM, N-CHANNEL FET	Q2
485527	MOSFET, IRF9523, .8 OHMS, 80V, P-CHANNEL	Q1
486716	IC, LM393N, DUAL COMPARATOR	IC3
515046	FUSE, 2 AMP, 250V, SLO-BLO	F301
515035	FUSE, 1/2 AMP, 250V, FAST-BLO	F302

NOTES:


- FOR TEST PROCEDURE SEE A2472-04.
- THE FOLLOWING COMPONENT IS NOT FITTED IN THIS ASSEMBLY: D6.

DO NOT SCALE UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES (mm) BREAK ALL SHARP EDGES TOLERANCES				TITLE		NOVAMATRIX MEDICAL SYSTEMS INC. WALLINGFORD, CT U.S.A. 06492			
				POWER SUPPLY BOARD ASSY, OXYPLETH - 520A					
DEC ± .082 (mm) ± .007 (mm) ± .003 (mm) ANG ±				DRAWN MBL 13/FEB/92		CHECKED AJE 4/2/92			
MATERIAL				MFG ENGR EH 5-15-92		APPROVED AEK 15-MAY-92			
FINISH				USED ON: D5873-01		SCALE: 2/1		SHEET 1 OF 1	
03	N403	9AUG96							
02	N341	10JAN94							
01	N316	31MAR93							
REV	R NO.	DATE							

This document contains information which is proprietary and the property of Novamatrix Medical Systems Inc. and may not be reproduced, used or disclosed in whole or in part without the written consent of Novamatrix Medical Systems Inc.



- △ NOTES:
1. UNLESS OTHERWISE SPECIFIED: RESISTANCE VALUES ARE IN OHMS, 1/4W, 1%.
 2. LAST REF DES USED: C7, D6, E302, F302, IC3, J302, L1, Q2, R16.
 3. CONNECT DGND AND ACND ONLY AT IC1-3.
 4. E301 AND D6 ARE FOR DC OPTION AND ARE NOT PART OF THIS ASSEMBLY.

DO NOT SCALE UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES (mm) BREAK ALL SHARP EDGES			SCHEMATIC, POWER SUPPLY BOARD, MODEL 515A SB, 520A, 1265 & 7100		 NOVAMATRIX Medical Systems, Inc. Wallingford, CT USA 06492			
MATERIAL:			DRAWN Tony Esposito Jul-10-91				CHECKED M.R. Lineback	
FINISH:			MFG. ENCR Ed Harvie 5-15-92		APPROVED A.E. Kersey 15-May-92			
03	N403	13Aug96	FILE: 24720303		SIZE	DWG NO.	CODE	REV.
02	N341	16Feb94	DATE: 13-Aug-96		D	2472	03	03
01	N316	1Apr93	TIME: 11:04:27		USED ON D2472-01		SCALE: NONE	
REV. "R" NO. BY DATE					SHEET 1 OF 1			

This document contains information which is proprietary and the property of Novamatrix Medical Systems Inc. and may not be reproduced, used or disclosed in whole or in part without the written consent of Novamatrix Medical Systems Inc.

REF DES	TYPE NO.	IC CHART														DEDCOUPLING CAP (.1uF)
		VCC	VCCI	VBACK	V IN	+VI	+V5	+V12	+VA	-V5	-V12	-VA	A GND	D GND	I GND	
IC1	DS232CPE	-	16	-	-	-	-	-	-	-	-	-	-	-	15	C2 TO I GND
IC2	78M05CKC	-	OUT	-	-	IN	-	-	-	-	-	-	-	-	GND	NONE
IC3	CNY17-III	-	-	-	-	-	-	-	-	-	-	-	-	-	-	NONE
IC4	CNY17-III	-	-	-	-	-	-	-	-	-	-	-	-	-	-	NONE
IC5	CNY17-III	-	-	-	-	-	-	-	-	-	-	-	-	-	-	NONE
IC6	CNY17-III	-	-	-	-	-	-	-	-	-	-	-	-	-	-	NONE
IC7	LM78L12ACZ	-	-	-	-	-	-	OUT	IN	-	-	-	GND	-	-	C13 TO A GND
IC8	LM79L12ACZ	-	-	-	-	-	-	-	-	-	OUT	IN	GND	-	-	C12 TO A GND
IC9	UC3843N	-	-	-	7 Δ	-	-	-	-	-	-	-	-	5	-	C16 TO D GND (.47uF)
IC10	CD4013B	-	-	14	-	-	-	-	-	-	-	-	-	7	-	C23 TO D GND
IC11	CD4093BE	-	-	14	-	-	-	-	-	-	-	-	-	7	-	C25 TO D GND
IC12	LP2950ACZ	-	-	OUT	IN	-	-	-	-	-	-	-	GND	-	-	C22 TO D GND (2.2uF)
IC13	MM74HC573N	20	-	-	-	-	-	-	-	-	-	-	-	10	-	C29 TO D GND
IC14	MM74HC573N	20	-	-	-	-	-	-	-	-	-	-	-	10	-	C30 TO D GND
IC15	MM74HC573N	20	-	-	-	-	-	-	-	-	-	-	-	10	-	NONE
IC16	HID641BOR1P6	33	-	-	-	-	-	-	-	-	-	-	-	1, 33	-	C33 TO D GND
IC17 Δ	HN27C101G-17	32	-	-	-	-	-	-	-	-	-	-	-	16	-	C40 TO D GND
IC18 Δ	PEEL18CVB	20	-	-	-	-	-	-	-	-	-	-	-	10	-	C49 TO D GND
IC19	MM74HC14N	14	-	-	-	-	-	-	-	-	-	-	-	7	-	C37 TO D GND
IC20	HM62256LP-12	-	-	28	-	-	-	-	-	-	-	-	-	14	-	C41 TO D GND
IC21	MM74HC14N	14	-	-	-	-	-	-	-	-	-	-	-	7	-	C48 TO D GND
IC22	MM74HC138	16	-	-	-	-	-	-	-	-	-	-	-	8	-	C44 TO D GND
IC23	MC74HC32N	14	-	-	-	-	-	-	-	-	-	-	-	7	-	NONE
IC24	MSM6242BRS	-	-	18	-	-	-	-	-	-	-	-	-	9	-	C42 TO D GND (2.2uF)
IC25	MC74HC32N	-	-	14	-	-	-	-	-	-	-	-	-	7	-	C43 TO D GND
IC26	TB820M	6	-	-	-	-	-	-	-	-	-	-	-	4	-	C52 TO D GND
IC27	SN76495AN	16	-	-	-	-	-	-	-	-	-	-	-	8	-	C56 TO D GND
IC28	MM74HC573N	20	-	-	-	-	-	-	-	-	-	-	-	10	-	C58 TO D GND
IC29	AD7528JN	17	-	-	-	-	-	-	-	-	-	-	1	5	-	C108 TO D GND
IC30	TL074CN	-	-	-	-	-	-	4	-	-	-	-	-	-	-	C83 TO A GND C84 TO A GND
IC31	TLC5491P	-	-	-	-	-	8	-	-	-	-	-	4	-	-	C70 TO A GND
IC32	LM317LZ	-	-	-	-	-	-	-	IN	-	-	-	-	-	-	NONE
IC33	MM74HC4051N	16	-	-	-	-	-	-	-	7	-	-	-	3	-	C68 TO D GND
IC34	CS5503JP	-	-	-	-	-	14 15 Δ	-	-	-	-	-	6 Δ 7	-	-	C75 TO D GND C74 TO D GND C76 TO D GND C77 TO D GND
IC35	LT1019CN8-2.5	-	-	-	-	-	-	2 Δ	-	-	-	-	4	-	-	C82 TO A GND
IC36	TL074CN	-	-	-	-	-	-	4	-	-	-	-	-	-	-	C71 TO A GND C86 TO A GND
IC37	CS5503JP	-	-	-	-	-	14 15 Δ	-	-	-	-	-	6 Δ 7	-	-	C79 TO D GND C78 TO D GND C80 TO D GND C81 TO D GND
IC38	AD712JN	-	-	-	-	-	-	8	-	-	-	-	-	-	-	C98 TO A GND C99 TO A GND
IC39	MM74HC4020N	16	-	-	-	-	-	-	-	-	-	-	-	8	-	C101 TO D GND
IC40	TLE2022CP	-	-	-	-	-	-	8	-	-	-	-	-	-	-	C91 C102
IC41	DG444DJ	12	-	-	-	-	-	13	-	-	-	-	5	-	-	C93 TO A GND C106 TO A GND
IC42 Δ	PEEL18CVB	20	-	-	-	-	-	-	-	-	-	-	-	10	-	C107 TO D GND
IC43	LM79L05ACZ	-	-	-	-	-	-	-	-	OUT	-	IN	GND	-	-	C65 TO A GND
IC44	LM78L05ACZ	-	-	-	-	-	-	OUT	-	IN	-	-	GND	-	-	C66 TO A GND
IC45	E1166	1	-	-	-	-	-	-	-	-	-	-	-	2	-	C110 TO A GND (100uF)
IC46	MC74HC32N	14	-	-	-	-	-	-	-	-	-	-	-	7	-	C47 TO D GND


NOTES:

- 1. IC17 IS PART OF S₀2 PROGRAM.
- 2. IC18 IS PART OF CSIO CONTROLLER PROGRAM.
- 3. IC42 IS PART OF TIMING SEQUENCER PROGRAM.
- 4. POWER IS INDIRECTLY CONNECTED THROUGH ADDITIONAL COMPONENTS.
- 5. UNLESS OTHERWISE SPECIFIED:
RESISTOR VALUES ARE IN OHMS, 1/4W, ±1%
- 6. LAST REFERENCE DESIGNATOR USED: C114, D18, F2, FB3, IC46, J109, L2, Q17, R64, RP13, T1, TP2, VR1, Y3
- 7. REFERENCE DESIGNATOR NOT USED: J103, LS1
- 8. TONE GENERATOR BOARD 2755-01 FITTED IN PLACE OF IC27.

		07 N403 6Aug96	DO NOT SCALE UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES (mm) BREAK ALL SHARP EDGES		TITLE	
		06 N362 23Mar95			SCHEMATIC, MAIN BOARD.	
		05 N348 5May94			MODEL 520A	
11 N681 12Aug99	04 N325 23Jul93	DEC ± (mm) _____ FRAC ± .007 - .003 (mm) (+.18 -.08) MATERIAL _____		DRAWN MRL 30APR92		
10 N557 3Jun98	03 N316 31MAR93			CHECKED AE 2APR92		
09 N557 28May98	02 N304 4SEPT92			MFC ENGR EH 29APR92		
08 N510 19Sep97	01 N302 2JUN92			APPROVED AK 29APR92		
REV R NO. BY DATE	REV R NO. DATE	FINISH _____		USED ON D2471-01		
				SIZE D	DRAWING NO. 2471	CODE 03
				SCALE: NONE	REV 11	
				SHEET 1 OF 5		



This document contains information which is proprietary and the property of Novamatrix Medical Systems Inc. and may not be reproduced, used or disclosed in whole or in part without the written consent of Novamatrix Medical Systems Inc.




To Sheets 2,3 and 4

11	N681	12Aug99	DO NOT SCALE UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES (mm) BREAK ALL SHARP EDGES	SCHEMATIC, MAIN BOARD, OXYPLETH - 520A	NOVAMATRIX Medical Systems, Inc. Wallingford, CT USA 06492
10	N557	3Jun98			
09	N557	27May98			
08	N510	26Sep97	MATERIAL:		
07	N403	13Aug96	FINISH:	DRAWN Tony Esposito Jan-30-1991	CHECKED M.R. Lineback 15May92
06	N362	5Apr95		MFG. ENGR Ed Harvie 5-15-92	APPROVED A.E. Kersey 15-May-92
05	N238	13Sep94	FILE: 24710311b DATE: 12-Aug-99 TIME: 15:40:55	USED ON D2471-01	
REV.	"R" NO.	BY DATE			
SIZE	DWG NO.	CODE	REV.		
D	2471	03	11		
SCALE:	NONE	SHEET	2 OF 5		

This document contains information which is proprietary and the property of Novamatrix Medical Systems Inc. and may not be reproduced, used or disclosed in whole or in part without the written consent of Novamatrix Medical Systems Inc.


CONTROL BUS

To Sheets 1, 3 and 4



11	N681	12Aug99	DO NOT SCALE UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES (mm) BREAK ALL SHARP EDGES
10	N557	3Jun98	
09	N557	27May98	
08	N510	26Sep97	MATERIAL:
07	N403	13Aug96	FINISH:
06	N362	5Apr95	
05	N348	13Sep94	FILE: 24710311c DATE: 17-Aug-99 TIME: 11:51:08
REV.	"R"	NO.	BY DATE


SCHEMATIC, MAIN BOARD, OXYPLETH - 520A	
DRAWN Tony Esposito Jan-30-92	CHECKED M.R. Lineback 15May92
MFG. ENGR Ed Harvie 5-15-92	APPROVED A.E. Kersey 15-May-92
USED ON	D2471-01

NOVAMATRIX 			
Medical Systems, Inc. Hallingford, CT USA 06492			
SIZE D	DWG NO. 2471	CODE 03	REV. 11
SCALE: NONE		SHEET 3 OF 5	

This document contains information which is proprietary and the property of Novamatrix Medical Systems Inc. and may not be reproduced, used or disclosed in whole or in part without the written consent of Novamatrix Medical Systems Inc.

DATA BUS
ADDRESS BUS
CONTROL BUS


To Sheets
1, 2 and 4




11	N681	12Aug99	DO NOT SCALE UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES (mm) BREAK ALL SHARP EDGES	SCHEMATIC, MAIN BOARD, OXYPLETH - 520A	NOVAMATRIX Medical Systems, Inc. Hartford, CT USA 06492	SIZE D	DWG NO. 2471	CODE 03	REV. 11
10	N557	3Jun98							
09	N557	27May98							
08	N510	26Sep97	MATERIAL:						
07	N403	13Aug96	FINISH:						
06	N362	5Apr95							
05	N348	13Sep94	FILE: 24710311d DATE: 12-Aug-99 TIME: 16:25:13	DRAWN Tony Esposito Jan-20-1991	CHECKED M.R. Lineback 22May92	APPROVED A.E. Kersey 22-May-92			
REV. "R" NO.	BY	DATE		USED ON D2471-01					

SCALE: NONE
SHEET 4 OF 5

This document contains information which is proprietary and the property of Novamatrix Medical Systems Inc. and may not be reproduced, used or disclosed in whole or in part without the written consent of Novamatrix Medical Systems Inc.



11	N681	12Aug99	DO NOT SCALE UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES (mm) BREAK ALL SHARP EDGES	SCHEMATIC, MAIN BOARD, OXYPLETH - 520A	NOVAMATRIX Medical Systems, Inc. Wallingford, CT USA 06492
10	N557	3Jun98			
09	N557	27May98		DRAWN Tony Esposito Jan-30-1991	CHECKED H.R. Lineback 22May92
08	N510	26Sep97	MATERIAL:	APPROVED A.E. Kersey 22-May-92	
07	N403	13Aug96	FINISH:	MFG. ENGR Ed Harvie 5-22-92	FILE: 24710311e DATE: 12-Aug-99 TIME: 17:40:59
06	N362	5Apr95		USED ON D2471-01	SIZE D
05	N348	13Sep94			DWG NO. 2471
REV. "R" NO.	BY DATE				CODE 03
					REV. 11
					SCALE: NONE
					SHEET 5 OF 5



SEE NOTE 1E

SECURE Y3 TO BOARD WITH FOAM TAPE 3/16 IN. X 1/16 IN.

SEE NOTE 1C

SEE NOTE 1C

SEE NOTE 1D

SEE NOTE 1C

2775-01: W1 ACROSS 520A & COM PINS
2775A-01: W1 ACROSS 2001 & COM PINS

NOTES:


1. GENERAL SPECIFICATIONS:
 - A. ASSEMBLY IS SHOWN FOR REFERENCE ONLY.
 - B. FOR TEST PROCEDURE SEE A2775-04.
 - C. WHEN MOUNTING IC2, Q5, AND Q16 ENSURE THAT THE HEIGHT IS NO HIGHER THAN C9.
 - D. HAND SOLDER TRANSFORMER T1. DO NOT WAVE SOLDER.
 - E. INSTALL 2755 BOARD AT IC27 WITH ORIENTATION AS SHOWN.

DRAWING TABULATION			
DWG NO.	REV	TITLE	NEXT ASSY
D2775-01	04 N820	MAIN BOARD ASSEMBLY, 520A	D5693-01
D2775A-01	04 N820	MAIN BOARD ASSEMBLY, 2001	D9400-01

This document contains information which is proprietary and the property of Novamatrix Medical Systems Inc. and may not be reproduced, used or in part without the written consent of Novamatrix Medical Systems Inc.

DO NOT SCALE UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES (mm)	DRAWN B. Lancaster DATE 28Jul99	
TOLERANCES .XX ±.01 [.254] .XXX ±.005 [.127] .XXXX ±.0025 [.064] FRAC ±1/64 ANG ±1/4°	CHECK MRL DATE 80c199	
04 N820 17Nov00	MFG ENGR T. Carroll DATE 80c199	TITLE MAIN BOARD ASSY
03 N804 11Oct00	ENGR DATE	SIZE D
02 N745 5Apr00	PROJ ENGR AJE DATE 80c199	DRAWING NO. 2775
01 N727 17Jan00	DATE 80c199	CODE 01
REV R NO. DATE	NEXT ASSY: SEE TAB	SCALE: 2/1
		SHEET 1 OF 2


SEE SHEET 2 FOR SOLDER SIDE ASSY




DRAWING TABULATION			
DWG NO.	REV	TITLE	NEXT ASSY
D2775-01	04 N820	MAIN BOARD ASSEMBLY, 520A	D5693-01
D2775A-01	04 N820	MAIN BOARD ASSEMBLY, 2001	D9400-01

This document contains information which is proprietary and the property of Novamatrix Medical Systems Inc. and may not be reproduced, used or disclosed in whole or in part without the written consent of Novamatrix Medical Systems Inc.

DO NOT SCALE UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES (mm)	DRAWN B. Lancaster DATE 29Jul99		
TOLERANCES .XX ±.01 [.254] .XXX ±.005 [.127] .XXXX ±.0025 [.064] FRAC 1/16" ANG ±1/4"	CHECK MRL DATE 8Oct99		
	MFG ENGR T. Carroll DATE 8Oct99	TITLE MAIN BOARD ASSY	
	ENGR DATE		
	PROJ ENGR AJE DATE 8Oct99	SIZE D	DRAWING NO. 2775
REV R NO.	DATE	CODE 01	REV SEE TAB
		SCALE: 2/1	SHEET 2 OF 2




SCHEMATIC, MAIN BOARD, 520A & 2001				Novamatrix			
Medical Systems, Inc. 1 Barnes Industrial Park Road Wallingford, CT USA 06492							
04	N820	17Nov00	Drawn: EJM	12Aug99	Checked: MRL	Size: B	Document Number: 2775 & 2775A
03	N804	5Oct00	Mfg: E7g	TC 07oct99	Approved: AJE 7Oct99	Code: 03	Rev: 04
02	N745	5Apr00	File: 27750304	Date: Friday, November 17, 2000		Sheet 1	of 4
01	N727	17Jan00					
Rev.	"R" No.	By Date					






- NOTES:
- UNLESS OTHERWISE SPECIFIED:
RESISTANCE VALUES ARE IN OHMS, 1/4 W, 1%
CAPACITANCE ARE IN FARADS
 - TONE GENERATOR BOARD 2755-01 FITTED IN PLACE OF IC27
 - LAST REFERENCE DESIGNATOR USED: C116, D18, F2, FB18, IC46, J109, JP2, L2, Q17, R68, RP13, SG2, T1, TP2, VR1, W1, Y3
 - REFERENCE DESIGNATORS NOT USED: C31, C36, D10, R23, R27
 - 2776-01 ASSY IS PART OF 2775A-01 ASSY.

SCHEMATIC, MAIN BOARD, 520A & 2001				Novamatrix Medical Systems, Inc. 1 Barnes Industrial Park Road Wallingford, CT USA 06492	
Drawn: EJM 12Aug99	Checked: MRL	Size: B	Document Number: 2775 & 2775A	Code: 03	Rev: 04
Mfg: N820 03	Eng: N804 02	File: TC 07oct99	Approved: AJE 7Oct99	Sheet 3	of 4
Rev.:	R* No.:	By Date:	Date: Friday, November 17, 2000		





NOTE: JUMPER W1 CONNECTS 2-3 FOR MODEL 520A, AND 1-2 FOR MODEL 2001

SCHEMATIC, MAIN BOARD, 520A & 2001				Novamatrix	
				Medical Systems, Inc.	
				1 Barnes Industrial Park Road	
				Wallingford, CT USA 06492	
04	N820	17Nov00	Drawn: EJM	12Aug99	Checked: MRL
03	N804	5Oct00	Eng: TC	07Oct99	Approved: AJE
02	N745	5Apr00	File: 27750304		
01	N727	17Jan00	Date: Friday, November 17, 2000		
Rev.	"R" No.	By Date			
			Size B	Document Number 2775 & 2775A	Code 03 Rev 04
					Sheet 4 of 4



NOTES:
1. COMPONENTS NOT FITTED: Q1, R11





			DO NOT SCALE UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES BREAK ALL SHARP EDGES		TITLE TONE GENERATOR REPLACEMENT BOARD ASSY		NOVAMETRIX  MEDICAL SYSTEMS INC. WALLINGFORD, CT U.S.A. 06492					
			.X ±.02 TOLERANCES .XX ±.01 .XXX ±.005 FRAC ±1/64 .XXXX ±.0010 ANG ±1/4°		DRAWN BL 5Jan98						CHECKED MRL 26Feb98	
			MATERIAL		MFG ENGR TJC 26Feb98		APPROVED RHD 26Feb98		SIZE B	DRAWING NO. 2755	CODE 01	REV 02
02	N699	40Oct99	~		USED ON: ~		SCALE: 4/1		SHEET 1 OF 1			
01	Prod	Release	FINISH									
REV	R NO.	DATE										




- NOTES:
- UNLESS OTHERWISE SPECIFIED, RESISTANCE VALUES ARE IN OHMS, 1/16W 1% CAPACITANCE VALUES ARE IN FARADS
 - LAST REF DES USED: C4, IC3, J2, Q1, R12
 - COMPONENTS NOT FITTED IN THIS ASSEMBLY: Q1, R11

			Schematic, Tone Generator Replacement Board		Novamatrix Medical Systems, Inc. 5 Technology Drive Wallingford, CT USA 06492		
			Drawn: R Daniels 23Dec97	Checked: MRL 25FEB98			
			Mfg Eng: TC 17FEB98	Approved: RHD 5FEB98	Size: B	Document Number: 2755	
02	N699	5Oct99	File: SCHEMATIC1			Code: 03	Rev: 02
01	Prod	Release	Date: Tuesday, October 05, 1999			Sheet: 1	of: 1
Rev.	"R" No	By Date					

This document contains information which is proprietary and the property of Novamatrix Medical Systems Inc. and may not be reproduced, used or disclosed in whole or in part without the written consent of Novamatrix Medical Systems Inc.



☆ TWIST BLACK & BROWN WIRES 6 TURNS CLOCKWISE.




NOTES:


1. BEFORE ASSEMBLING ITEMS 1 & 4, APPLY 1 DROP OF ITEM 2 TO DRESS BEZEL IN AREA B AS SHOWN IN DETAIL A AND TO AREA D IN REAR VIEW C-C. FINGER TIGHTEN NUT ONTO CONNECTOR UNTIL IT BOTTOMS OUT AGAINST ITEM 1. ALLOW ASSY TO CURE FOR 24 HOURS.
2. WIRE WILL BE PREPPED WHEN CONTACT (ITEM 6) IS INSTALLED. FOR REFERENCE ONLY: 1/8 STRIP DO NOT TIN.
3. FOR EPOXY MIXING ASSY PROCEDURE SEE A9779-33.

			DO NOT SCALE UNLESS OTHERWISE SPECIFIED. DIMENSIONS ARE IN INCHES (mm). BREAK ALL SHARP EDGES.		TITLE		NOVAMETRIX MEDICAL SYSTEMS INC. WALLINGFORD, CT U.S.A. 06492				
			TOLERANCES		CABLE ASSY, S _o O ₂ INPUT, MODEL 520A						
05	N718	25Jan00	DEC ±	(mm)	DRAWN	MRL	CHECKED	SIZE	DRAWING NO.	CODE	REV
04	N681	16Jul99	FRAC ±	1/16 (mm)	23.OCT.91	A.J.E.		C	5728	01	05
03	N366	17May95	HOLES ±	+0.07 -0.03 (mm) (+.18 -.08)	MFG ENGR	B.N.	APPROVED				
02	N302	26MAY92	ANG ±		12FEB92	2FEB92					
01	N299	26MAR92	FINISH								
REV	R NO.	DATE				USED ON: D5693-01		SCALE: 2/1		SHEET 1 OF 1	


This document contains information which is proprietary and the property of Novamatrix Medical Systems Inc. and may not be reproduced, used or disclosed in whole or in part without the written consent of Novamatrix Medical Systems Inc.



- △ NOTES:
- UNLESS OTHERWISE SPECIFIED: RESISTANCE VALUES ARE IN OHMS, 1/4W, 1%.
 - LAST REFERENCE DESIGNATOR USED: D4, E2, P105, R2.

			DO NOT SCALE UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES (mm) BREAK ALL SHARP EDGES	SCHEMATIC, ALERT BOARD, MODEL 520A		NOVAMATRIX  Medical Systems, Inc. Wallingford, CT USA 06492			
			MATERIAL: _____						DRAWN Tony Esposito 1-May-92
			FINISH: _____	MFG. ENGR Ed Harvie 5-15-92	APPROVED A.E. Kersey 15-May-92				
01	N304	9Nov92	FILE: 24730301 DATE: 18-Apr-96 TIME: 16:30:05	USED ON D2473-01		SCALE: NONE		SHEET 1 OF 1	
REV.	"R" NO.	BY DATE							

This document contains information which is proprietary and the property of Novometrix Medical Systems Inc. and may not be reproduced, used or disclosed in whole or in part without the written consent of Novometrix Medical Systems Inc.








SCHMATIC DIAGRAM

★ S11 & S12 ARE ALLOWED FOR IN THE CIRCUITRY BUT ARE NOT USED FOR THIS KEYPANEL.

DO NOT SCALE UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES (mm) BREAK ALL SHARP EDGES			TOLERANCES		TITLE	
05	N388	2Apr98	DEC ±	(mm)	MEMBRANE KEYPANEL, OXYPLETH - 520A	
04	N316	25MAR93	FRACTION ±	(mm)	NOVOMETRIX MEDICAL SYSTEMS INC. WALLINGFORD, CT U.S.A 06492	
03	N296	3/31/92	HOLES ±	(mm) (+.18 - .08)	DRAWN MRL SOMAY91	CHECKED A.P.
02	N296	3/11/92	ANG ±	(mm)	MFG ENGR	APPROVED A.P.
01	N292	2/19/92	FINISH			
REV	R NO.	DATE	USED ON: D5719-01		SCALE: NONE	SHEET 2 OF 2

This document contains information which is proprietary and the property of Novamatrix Medical Systems Inc. and may not be reproduced, used or disclosed in whole or in part without the written consent of Novamatrix Medical Systems Inc.



- NOTES:
- 1 APPLY A SMALL DROP OF ADHESIVE, ITEM 23, TO STUDS, ITEM 13, BEFORE ASSEMBLY.
 - 2 APPLY A SMALL DROP OF ADHESIVE, ITEM 23, TO SCREWS, ITEM 18, BEFORE ASSEMBLY.


DO NOT SCALE UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES (mm) BREAK ALL SHARP EDGES				TITLE		NOVAMATRIX			
TOLERANCES				REAR PANEL ASSY, MODEL 520A AND 515A SB		MEDICAL SYSTEMS INC. WALLINGFORD, CT U.S.A. 06492			
DEC ±	FRAC ±	HOLE ±	MFG ±	DRAWN	CHECKED	SIZE	DRAWING NO.	CODE	REV
03	N325	BJU193	---	MRL	AE	D	5673	01	03
02	N316	24MAR93	---	EH	AK				
01	N304	14OCT92	FINISH	---	---				
REV	R NO.	DATE	FINISH	---	---				

USE ON: 5693-01/5750-01





SCALE: 1/1

SHEET 1 OF 1

This document contains information which is proprietary and the property of Novamatrix Medical Systems Inc. and may not be reproduced, used or disclosed in whole or in part without the written consent of Novamatrix Medical Systems Inc.



REAR VIEW A-A





WIRING DIAGRAM
NOTE: 1/8 STRIP & TIN ON EACH END FOR ALL WIRES.

NOTES:

- ① COVER WIRE CONNECTIONS WITH HEAT SHRINK, ITEM 13, 5 PL @ 3/8 L.
- ② COVER LEAD OF RESISTOR, ITEM 4 FROM P1-25 WITH 3/8 L TUBING, ITEM 11.

DO NOT SCALE UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES (mm) BREAK ALL SHARP EDGES				TITLE RS-232-C TEST FIXTURE ASSY, MODEL 515 & 515A			
TOLERANCES DEC ± ~ (mm) (± ~) FRAC ± 1/32 (mm) (± ~) HOLES +.007 -.003 (mm) (+.18 -.08) ANG ± ~				DRAWN RGG 4Dec95			
MATERIAL ~				MFG ENGR EBE 3-20-96		APPROVED ASB 18-Mar-96	
01	N338	4Dec95	FINISH	~	USED ON: ~	SCALE: 2/1	SHEET 1 OF 1




Mark the test jack part number, as shown, on the completed assembly

Parts List

Qty	Part#	Description
1	152014	Capacitor, 1uF
1	211705	Connector, 7 Pin, Hypertronics
1	472026	Resistor, 6.81K
1	472030	Resistor, 10K
1	472058	Resistor, 100K
1	481000	Diode, Zener, ICL8069
4	481501	Diode, 1N4148

Connector Pin Out
Rear View



Novamatrix Medical Systems		
Title Front End Saturation Test Jack		
Size A	Document Number 5942-00	Rev 00
Date:	Wednesday, June 28, 2000	Sheet 1 of 1